

TRENDELENBURG POSITIONING SYSTEM (TPS)

Instructions For Use

INTENDED USE

The Trendelenburg Positioning System is used to help prevent sliding of patients on the surgical table placed in the Trendelenburg position while assisting in protection from nerve damage and pressure ulcers.



Figure 1:
Aligning pad on table




Figure 2:
Securing straps to
accessory rail




Figure 3:
Placing lift sheet
across pad




Figure 4:
Patient position




Figure 5:
Securing patient's
arms with lift sheets





Figure 6:
Attaching body strap
to accessory rail



Figure 7:
Attaching the body
straps across the
patient's chest



Step 6. Position the patient on the pad

Place the patient onto the pad in a position that will allow tucking of the arms next to the patient's body. Make sure that the shoulders do not extend past the edge of the pad and that the skin at the scapula and sacrum are in direct contact with the surface of the pad (refer to Figure 4).

Step 7. Secure the patient's arms with the lift sheet

Wrap the lift sheet around the patient's arm. Tuck the end of the lift sheet under the patient (refer to Figure 5). Do not tuck the arms too tightly. Be sure to check the position of the patient's fingers around the surgical table mechanisms and for proper circulation.

Step 8. Secure the patient with the body strap

- Wrap the lift sheet around the patient's arm. Tuck the end of the lift sheet under the patient (refer to Figure 5). Do not tuck the arms too tightly. Be sure to check the position of the patient's fingers around the surgical table mechanisms and for proper circulation.
- Move to the opposite side of the surgical table. Locate the body strap with the large rectangle of hook-and-loop material. Repeat the above step, holding the body strap with the large rectangle of hook-and-loop material facing the patient.
- Affix the two (2) body straps to each other so that they are secure across the patient's chest (refer to Figure 7). Do not over tighten.

MAINTENANCE

Inspect the pad and straps. Do not use products with damage or deterioration.

STORAGE AND TRANSPORT

The packaged products shall be stored in a cool and dry location. Keep away from direct sunlight, corrosive substances and in a well-ventilated indoor location.

The pads should always be stored and handled flat. Do not store the pads rolled or folded.



Not made with
natural rubber latex



EUROPEAN REPRESENTATIVE
GRI-Alleset Limited B.V.
Kingsfordweg 151, 1043 GR
Amsterdam, THE NETHERLANDS



MANUFACTURER
GRI Medical & Electronic Technology Co., Ltd.
www.alleset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA



IMPORTER
GRI-Alleset Limited B.V.
Kingsfordweg 151, 1043 GR
Amsterdam, THE NETHERLANDS

U.S. HEADQUARTERS
Alleset, Inc.
www.alleset.com
4142 Industry Way
Flowery Branch, GA 30542 U.S.A.