VERESS NEEDLE

Instructions for Use

VERESS NEEDLE

RX ONLY

SINGLE PATIENT USE; DO NOT REUSE

NOT MADE WITH NATURAL RUBBER LATEX

PHTHALATE FREE

IMPORTANT INFORMATION

Please read all instructions, warnings, and precautions before use. Correct application is essential for proper functioning of the product. Use only by a physician or properly licensed practitioner who is familiar with the purpose for which it is intended and only for the use it was intended.

INTENDED USE

The Veress Needle is intended to be used in minimally invasive surgical procedures in conjunction with a pump to allow for inflation of the pneumoperitoneum. This device will be used according to the physician's decision to utilize a minimally invasive procedure and may be used on adult patients of any ethnicity and where the health of the patient will allow for a minimally invasive surgical procedure as determined by the physician.

CONTRAINDICATIONS

This instrument is not intended for use when minimally invasive techniques are contraindicated.

PRECAUTIONS

This device is designed for single use only. Do not reuse, reprocess, or re-sterilize. Do not use if the package is opened or damaged.

- Minimally invasive procedures should be performed only by physicians having adequate training and knowledge of minimally invasive techniques and according to the medical facility's procedure(s).
- Verify compatibility with all instrumentation, devices and accessories from different manufacturers prior to the initiation of the procedure.
- Do not attempt to use the Veress Needle if the "red safety indicator" in the handle does not slide in the direction of the stopcock.
- After removing the Veress Needle, always inspect the site for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device and other instruments.

WARNINGS

Entry into the abdominal cavity using a sharp instrument can cause major vessel injury and damage to adjacent vital intraperitoneal structures. Before beginning any endoscopic procedure, the surgeon must be familiar with the anatomy pertinent to the procedure and perhaps unique to the individual patient. Although the blunt stylet safety mechanism is intended to reduce the risk of injury related to penetrating vital intraperitoneal structures, correct technique with respect to insertion of the Veress Needle is essential. To further minimize the risk of injury be sure to:

- Locate important anatomical objects such as the aorta and both iliac crests to properly place the Veress Needle.
- Place the needle with a careful smooth motion in order to avoid impaling structures against the sacral promontory.
- Maintain awareness of the patient's position on the operating table

This device is intended for single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may result in the following complications:

- Compromised structural integrity and/or device failure that in turn may result in patent injury, illness or death.
- Risk of contamination and/or patient infection or crossinfection that in turn may lead to injury, illness, or death of the patient. Maintain awareness of the patient's position on the operating table.
- A maximum intra-abdominal pressure (IAP) of 30mmHg is recommended for use with this Veress Needle.

Products delivered non-sterile which are ultimately for use in sterile environments are intended to be sterilized by the sterile packager. This includes further packaging and sterilization according to the sterile packagers validated processes.

In the event of a patient injury related to the device please report the details to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

STERILIZATION AND REUSE INFORMATION

These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by GRI-Alleset. Other methods of sterilization, multiple sterilizations, or reuse could cause a device malfunction and possibly result in injury to the patient.

INSTRUCTIONS FOR USE

Important: This package insert is designed to provide directions for use of the Veress Needle. It is not a reference to minimally invasive techniques.

- Using sterile technique, remove the instrument from the package. Remove the cap from the end of the needle and dispose of it. To avoid damage, do not flip the instrument into the sterile field.
- Make a small incision to insert the Veress Needle. Hold the Veress Needle handle between the thumb and index finger with the middle finger supporting the rest of the handle.
- During insertion, inspect the instrument handle to ensure that the "red safety indicator" fully exposes itself beyond the plastic ergonomic handle. This action retracts the blunt stylet and exposes the sharp needle tip to allow for penetration. Do not attempt to use the Veress Needle if the blunt stylet does not retract into the needle sleeve. Once the blunt stylet is free of tension from the tissue, the "red safety indicator" resets itself back into the distal portion of the handle.
- Use appropriate technique to ensure that the Veress Needle is safely in the abdominal cavity. For example, if techniques such as the manometer test or the saline drop test indicates correct placement, insufflate the abdomen after attaching the Veress tubing/filter set. A maximum intra-abdominal pressure (IAP) of 30mmHg is recommended for use with this Veress Needle.
- Once the appropriate pneumoperitoneum pressure has been reached, remove the Veress Needle from the abdomen and dispose of the Veress Needle per the facility's sharps disposal protocol.

MAINTENANCE

During use keep the lumen of the Veress Needle clean and free of air flow obstructions. Inspect the Veress Needle before and during use, and do not use products with damage or deterioration.

SHELF LIFE

The shelf life of the product is five (5) year if stored under general conditions.

STORAGE AND TRANSPORT


The packaged products shall be stored in a cool, dry location.

DISPOSAL

Dispose of the product as sharp/biohazard waste in accordance with local municipal, provincial, or hospital policies.

ILLUSTRATION AND NOMENCLATURE

- Needle
- Luer Lock Connector
- Valve
- Handle
- Red Safety Indicator
- Blunt Stylet

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GRI Medical & Electronic Technology Co., Ltd.
Emergo Europe B.V.
Alleset Inc.
U.S.A.