Instructions For Use INSUFFLATION TUBING

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
Insufflationsslang

Part No. 37-1000-S, 37-1003-S, 37-2000-S, 37-2003-S, 27-2005-S, 37-2022-S

EC REP EUROPEAN REPRESENTATIVE
Emergo Europe B.V.
Westervoortsedijk 60, 6827, AT Arnhem, The Netherlands

IMPORTER
GRI-Alleset Limited B.V.
Kingsfordweg 151, 1043 GR, Amsterdam, THE NETHERLANDS

U.S. HEADQUARTERS
Alleset, Inc.
www.alleaset.com
4142 Industry Way, Flowery Branch, GA 30542 U.S.A.

MANUFACTURER
GRI Medical & Electrical Technology Co., Ltd.
www.alleaset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

SM-IFU-InsufflationTubing-S V1 20240710

EN

Insufflation Tubing Instructions For Use

IMPORTANT INFORMATION

Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of product.

INTENDED USE

Insufflation tubing are intended to be used by a physician or technician in certain minimally invasive surgical procedures to assist with transferring the gas needed to establish a pneumoperitoneum in the patient. The product is intended to be used sterile and is single-use. Patient use will be determined by the attending physician.

INTENDED USERS

Physician and technician

USE ENVIRONMENT

OR surgical environment and are to be used sterile and single use.

INTENDED PATIENT POPULATION


These products may be used on patient's that will require minimally invasive surgical procedures and are not restricted by age, gender or ethnicity except as determined by the attending physician.

INDICATIONS

These products may be used when the physician determines that minimally invasive surgery is required.

WARNINGS

- Insufflation tubing is not intended to be used for introduction of fluid into the body.
- This system is only recommended for laparoscopic surgery
- Tubing sets will be connected to a trocar or veress needle and should not be used or inserted directly into the body, organs, and tissues.
- Inspect the product and packaging for visible cracks or broken components.
- DO NOT use if filter or housing is damaged.
- Purge device of air prior to delivering gas to patient.
- In the event of a patient injury related to the device please report the details to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

EC REP EUROPEAN REPRESENTATIVE
Emergo Europe B.V.
Westervoortsedijk 60, 6827, AT Arnhem, The Netherlands

IMPORTER
GRI-Alleset Limited B.V.
Kingsfordweg 151, 1043 GR, Amsterdam, THE NETHERLANDS

MANUFACTURER
GRI Medical & Electrical Technology Co., Ltd.
www.alleaset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

U.S. HEADQUARTERS
Alleset, Inc.
www.alleaset.com
4142 Industry Way, Flowery Branch, GA 30542 U.S.A.

CAUTIONS

DO NOT reuse or attempt to re-sterilize. Usage of this device other than directed by a physician or in these instructions may result in device malfunction and/or patient injury.

DIRECTIONS

1. Open the insufflation tubing set in the sterile field using aseptic technique.
2. The adapter or luer lock closest to the filter should be connected to the insufflator.
3. This connection should be made with the luer lock or an adapter that is compatible with the insufflator.
4. The opposite end of the tubing should be connected to either a trocar or veress needle (pneumo needle) for initial set-up. This connection should be made with the luer lock or an appropriate adapter.
5. Immediately before delivering gas to the patient, purge air from tubing sets and instruments with CO2.
6. If using a veress needle, move the insufflation tubing from the veress needle to the fitting on the trocar after pneumoperitoneum is established.

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
MANUFACTURER

GRI Medical & Electrical Technology Co., Ltd.
www.alleset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

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
MANUFACTURER

GRI Medical & Electrical Technology Co., Ltd.

www.alleset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

EC

REP

EUROPEAN REPRESENTATIVE

Emergo Europe B.V.
Westervoortsedijk 60, 6827, AT Arnhem, The Netherlands

IMPORTER

GRI-AlleSet Limited B.V.
Kingsfordweg 151, 1043 GR, Amsterdam, THE NETHERLANDS

U.S. HEADQUARTERS

AlleSet, Inc.
www.alleset.com
4142 Industry Way, Flowery Branch, GA 30542 U.S.A.

MANUFACTURER

GRI Medical & Electrical Technology Co., Ltd.

www.alleset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

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
GRI-Alleaset Limited B.V.
Kingsfordweg 151, 1043 GR, Amsterdam, THE NETHERLANDS

U.S. HEADQUARTERS

Alleaset, Inc.
www.alleaset.com
4142 Industry Way, Flowery Branch, GA 30542 U.S.A.

MANUFACTURER

GRI Medical & Electrical Technology Co., Ltd.

www.alleaset.com
1805 Honggao Road, 314031
Jiaxing, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

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