

# gentleheel<sup>®</sup>

Gently Precise Heel Incisions  
For Newborns & Toddlers



**gentleheel<sup>®</sup>** devices offer a safe and efficient way for acquiring blood samples from the heels of newborns and toddlers. Utilizing a sweeping incision, **gentleheel<sup>®</sup>** devices are engineered with features that help ensure ease of use and minimal discomfort.

## Intuitive Features

- Incision technique is less traumatic than puncture
- Permanent blade retraction minimizes possible reuse or injury
- Welded plastic casing prevents accidental exposure to blade
- Color coded for easier size selection
- Hardened stainless steel blade improves incision accuracy
- Multiple visual indicators reinforce incision site
- Beveled blade generates a smoother, less traumatic incision
- Ergonomic design for simplified operation



CODE	PKG	SIZE	COLOR	DESCRIPTION
TGH10X50	50 EA/PK (10 PK/CS)	2.00 mm / 3.00 mm	Blue	gentleheel <sup>®</sup> , Toddler
NGH10X50	50 EA/PK (10 PK/CS)	1.00 mm / 2.50 mm	Green	gentleheel <sup>®</sup> , Newborn
PGH10X50	50 EA/PK (10 PK/CS)	.85 mm / 1.75 mm	Lavender	gentleheel <sup>®</sup> , Premie
MPGH10X50	50 EA/PK (10 PK/CS)	.65 mm / 1.40 mm	Yellow	gentleheel <sup>®</sup> , Micro-Premie