

# Protective Apparel

## ComfortGuard® Multilayered SMS Light Weight Isolation Gown

	CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
	60-3040	Multilayered SMS Light Weight Isolation Gown	L	100	21x14x18" (54x35x46cm)
	60-3041	Multilayered SMS Light Weight Isolation Gown	XL	100	21x14x18" (54x35x46cm)

### Materials

- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency

### Product Attributes


- Full coverage
- Regular and extra-long lengths
- Tab neck
- Tie waist
- Elastic cuffs
- Sonic bonded seam construction

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.



Code	A: Total Length	B: Chest Length	C: Sleeve Length
60-3040	46" (116.8cm)	28.5" (72.3cm)	22.5" (57.1cm)
60-3041	48" (121.9cm)	30.5" (77.4cm)	23.5" (59.6cm)
Tolerance	±0.5" (1.3cm)	±1" (2.54cm)	±1" (2.54cm)

### Product Certifications/Standards

- ANSI/AAMI PB70:2022 Level 1 Liquid Barrier Protection
- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/ Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

