

# CORE Room Turnaround Kits<sup>®</sup>

## Kick Bucket Bags

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
50-KBBC01	Kick Bucket Bag, Clear, 1.5mil, 16x24" (41x61cm)	300	12x10x9" (31x25x22cm)
50-KBBR01	Kick Bucket Bag, 0.0013", Red w/ black printing, 23x24" (58x61cm)	250	18x12x7" (46x30x18cm)
50-KBBR02	Kick Bucket Bag, 0.0013", Red printed, 24x31" (61x79cm)	350	18x12x7" (46x30x18cm)
50-KBBB01	Kick Bucket Bag, 0.0013", Red w/ black printing, 24x31" (61x79cm)	100	16x12.5x10" (40x32x25cm)
KBRD2431	Kick Bucket Bag, 0.0013", Red w/ black printing, 24x31" (61x79cm)	350	18x12x7" (46x30x18cm)
KBRD2324	Kick Bucket Bag, 0.0013", Red w/ black printing, 23x24" (58x61cm)	500	18x12x7" (46x30x18cm)

### Materials

- Clear Polyethylene (PE) Film
- Red Polyethylene (PE) Film

### Product Attributes

- Red Kick Bucket Bag Includes Printed Bio-hazard Warning

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations


- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



50-KBBC01



Alleset  
4142 Industry Way  
Flower Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-002 2022  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



GRI-Alleset Limited B.V.  
Kingsfordweg 151,  
1043 GR Amsterdam,  
The Netherlands  
www.alleset.com

