

Surgical Accessories

Cautery Tip Cleaners

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
25-0200	Cautery Tip Cleaner, 2x2x0.2" (5x5x0.6cm)	2000	28x18x7" (65x45x17cm)
25-0200-S	Cautery Tip Cleaner, 2x2x0.2" (5x5x0.6cm)	100	11x9x5" (29x23x13cm)

Materials

- Foam Pad
- Adhesive Tape
- X-ray Detectable Element


Product Attributes

- Adhesive backing for universal placement
- X-ray detectable element for added safety

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.



Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards

- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex IX, Rule 1
- Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

