

Facial Protection

Eye Shields

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
33-0500	Eye Shields, Pre-Assembled, Clear w/ Multi-colored Bands	200	18x14x15" (46x35x38cm)
33-0501	Eye Shields, Lens, Clear w/ Multi-colored Bands	2500	17x10x9" (43x25x22cm)
33-0502	Eye Shields, Frames, Polypropylene, Black	200	17x11x16 (43x29x40cm)

Materials

- Frames: Polypropylene (PP)
- Lenses: PET (Polyethylene Terephthalate)

Product Attributes

- Fits over prescription glasses
- Reusable frame
- Removable protective film to protect lens until point of use
- Anti-fog lens

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 2 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 2 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards



- CE certified as Category 1 PPE per Regulation (EU) 2016/42

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

