

Facial Protection

Face Shield

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
 33-1001	Face Shield, Full Length	50	15x13x11" (37x32x27cm)

Materials

- Lens: PET (Polyethylene Terephthalate)
- Protective Film- Clear Polyethylene (PE)
- Elastic Strap

Product Attributes

- Full coverage with protective film liner
- Anti stat coating for fog free vision

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Product Certifications/Standards


- Class IIb, Sterile, according to MDD93/42/EEC Annex IX, Rule

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

