

Basins and Trays

Graduated Pitchers

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
15-1200	Pitcher, Graduated, 1200cc	70	19x13x16" (49x34x41cm)
15-1400	Pitcher, Graduated, 1400cc	160	26x26x21" (65x65x54cm)

Materials

- Polypropylene

Product Attributes

- Pour spout for easy transfer of liquids
- Integral handle design offers stability

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 2

Manufacturing Certifications/Registrations


- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



15-1200




15-1400



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

