

# ComfortGuard®

## Multilayered SMS Medium Weight Isolation Gown

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-3050	Multilayered SMS Medium Weight Isolation Gown	L	100	21x13x14" (53x35x36cm)
60-3051	Multilayered SMS Medium Weight Isolation Gown	XL	100	21x13x14" (53x35x36cm)

### Materials

- Multilayered SMS Material with Antistatic Properties

### Product Attributes

- Adjustable tape tab closure
- AAMI Level II Protection Level

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- ANSI/AAMI PB70:2012 Level 2 Liquid Barrier Protection
- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Code	A: Arm Hole Width	B: Total Length	C: Chest Length	D: Tie Length
60-3050	10.5" (26.6cm)	46" (116.8cm)	28.5" (72.3cm)	68" (172.7cm)
60-3051	11.5" (29.2cm)	48" (121.9cm)	30.5" (77.4cm)	70" (177.8cm)
Tolerance	±0.5" (1.3cm)	±1" (2.54cm)	±1" (2.54cm)	±1" (2.54cm)



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

