

# Surgical Sharps Safety

## Magnetic Instruments Pad

	CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
	31-1016	Instrument Pad, Magnetic, 10x16" (25.4x40.6cm)	100	8x17x21" (20x44x54cm)
	31-1016-S	Instrument Pad, Magnetic, 10x16" (25.4x40.6cm)	30	15x13x7" (38x33x18cm)
	31-2016	Instrument Pad, Magnetic, 16x20" (40.6x50.8cm)	50	8x17x21" (20x44x54cm)
	31-2016-S	Instrument Pad, Magnetic, 16x20" (40.6x50.8cm)	20	15x13x7" (38x33x18cm)

### Materials

- Blue Polyurethane (PU) Foam
- Magnet (Iron and Polyethylene)
- Blue Polyethylene (PE) Film

### Product Attributes

- Magnet securely holds multiple instruments

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
[www.alleset.com](http://www.alleset.com)  
 SM-TDS-002 2022  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
 Westervoortsedijk 60,  
 6827 AT Arnhem,  
 The Netherlands  
[www.emergobyul.com](http://www.emergobyul.com)

