

Patient Gowns

Blue

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
60-7001	Patient Robe, Blue	100	16.5x11x12.5" (42x28x32cm)
60-7002	Patient Robe, Blue	100	20x14x21" (52x35x54cm)

Materials

- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency

Product Attributes

- Full coverage provides modesty

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



60-7001



60-7002

Code	A: Back Length	B: Sleeve Length
60-7001	45" (114cm)	10" (25cm)
60-7002	49" (125.5cm)	14" (35cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

