

Protective Apparel

Bouffant Caps

CODE	DESCRIPTION (AxB)	CASE QTY	CASE DIMENSIONS
60-1000	Bouffant Cap, Blue, 21" (53cm)	1000	25x23x13" (64x59x32cm)
60-1001	Bouffant Cap, Blue, 24" (61cm)	750	25x23x13" (64x59x32cm)
60-1010	Bouffant Cap, White, 21" (53cm)	300	16x10x12" (41x25x30cm)
60-1011	Bouffant Cap, White, 24" (61cm)	300	16x10x12" (41x25x30cm)
60-1030	Protective Bouffant Cap, Blue, 24" (61cm)	700	27x21x12" (70x54x30cm)

Materials

- Polypropylene (PP)
- Heavy-Duty Spunbond Polypropylene (PP)
- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency (60-1030)

Product Attributes

- Lightweight and breathable for dry particle protection (60-1000, 1001, 1010, 1011)
- Tri-layer material for light/moderate protection (60-1030)

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards



- CE certified as Class 1, according to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified.

Storage

- General warehouse conditions



CODE	A: Length/Width of Bouffant	B: Length/Width of Elastic
60-1000	12" (30cm)	5" (12.7cm)
60-1001	13" (33cm)	5" (12.7cm)
60-1010	12" (30cm)	5" (12.7cm)
60-1011	13" (33cm)	5" (12.7cm)
60-1030	13" (33cm)	5" (12.7cm)
Tolerance	±1" (2.54cm)	±1" (2.54cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-0032023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.alleset.com