

# Protective Apparel

## Shoe Covers

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-2001	Regular Duty Shoe Cover	Universal	400	22x13x11" (56x34x28cm)
60-2002	Regular Duty Shoe Cover	XL	400	22x13x11" (56x34x28cm)
60-2003	Regular Duty Shoe Cover	Universal	500	28x13x11" (70x34x28cm)
60-2010	Regular Duty Shoe Cover w/ Traction	Universal	400	22x13x11" (56x34x28cm)
60-2011	Regular Duty Shoe Cover w/ Traction	XL	400	22x13x11" (56x34x28cm)
60-2020	Heavy Duty Shoe Cover	Universal	400	22x13x11" (56x34x28cm)
60-2030	Heavy Duty Shoe Cover w/ Traction	Universal	400	22x13x11" (56x34x28cm)
60-2031	Heavy Duty Shoe Cover w/ Traction	Universal	500	28x13x11" (70x34x28cm)
60-2032	Heavy Duty Shoe Cover w/ Traction	XL	400	22x13x11" (56x34x28cm)
60-2033	Heavy Duty Shoe Cover w/ Traction	XL	500	28x13x11" (70x34x28cm)
60-2034	Heavy Duty Shoe Cover w/ Traction	2XL	400	22x13x11" (56x34x28cm)
60-2035	Heavy Duty Shoe Cover w/ Traction	2XL	500	28x13x11" (70x34x28cm)

### Materials

- Regular Duty Shoe Cover: Polypropylene (PP)
- Heavy-Duty Shoe Cover: Multilayered SMS Material with Antistatic Properties

### Product Attributes

- Available with and without traction
- Dispenser box or bulk packing options

### Sterilization

**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as Category 1 PPE per Regulation (EU) 2016/425

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



CODE	A: Length	B: Back Height
60-2001	15.5" (39cm)	6.5" (16.5cm)
60-2002	16.5" (42cm)	6.5" (16.5cm)
60-2003	15.5" (39cm)	6.5" (16.5cm)
60-2010	15.5" (39cm)	6.5" (16.5cm)
60-2011	16.5" (42cm)	6.5" (16.5cm)
60-2020	15.5" (39cm)	6.5" (16.5cm)
60-2030	15.5" (39cm)	6.5" (16.5cm)
60-2031	15.5" (39cm)	6.5" (16.5cm)
60-2032	16.5" (42cm)	6.5" (16.5cm)
60-2033	16.5" (42cm)	6.5" (16.5cm)
60-2034	17.5" (44cm)	7" (18cm)
60-2035	17.5" (44cm)	7" (18cm)
Tolerance	±0.8" (2cm)	±0.5" (1.3cm)



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-002 2022  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
www.emergobyul.com

