

Protective Apparel

Apron

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
60-8000	Apron, Blue, 57x42" (146x106cm)	80	17x11x6" (42x30x17cm)

Materials

- Blue Polyethylene (PE) Film

Product Attributes

- Heavyweight film for high fluid procedures

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Establishes at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Category 1 PPE per Regulation (EU) 2016/425

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Item	A	B	C	D
60-8000	57" (145cm)	42" (107cm)	37" (94cm)	20" (51cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Prinsessegracht 20
2514 AP The Hague
The Netherlands
www.emergobyul.com

