

# Surgical Sharps Safety

## SECURE Needle Counters

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
<b>N</b> 33-01037	10-Count, Foam Block, Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	15x13x7" (38x33x18cm)
<b>S</b> 33-01037-S	10-Count, Foam Block, Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	15x13x7" (38x33x18cm)
<b>N</b> 33-01037Y	10-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	15x13x7" (38x33x18cm)
<b>S</b> 33-01037Y-S	10-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	15x13x7" (38x33x18cm)
<b>N</b> 33-01580	15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x5" (5x13cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-01580-S	15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x5" (5x13cm)	100	13x10x15" (33x25x38cm)
<b>N</b> 33-01580Y	15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x5" (5x13cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-01580Y-S	15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow, 2x5" (5x13cm)	100	13x10x15" (33x25x38cm)
<b>N</b> 33-01585	15-Count, Foam Strips, Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-01585-S	15-Count, Foam Strips, Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	13x10x15" (33x25x38cm)
<b>N</b> 33-01585Y	15-Count, Foam Strips, Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-01585Y-S	15-Count, Foam Strips, Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	13x10x15" (33x25x38cm)
<b>N</b> 33-02005	20-Count, Foam Block, Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	14x9x20" (36x22x50cm)
<b>S</b> 33-02005-S	20-Count, Foam Block, Magnet, Blade Remover, Red, 2x5" (5x13cm)	80	13x15x15" (33x38x38cm)
<b>N</b> 33-02005Y	20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm)	100	14x9x20" (36x22x50cm)
<b>S</b> 33-02005Y-S	20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm)	80	13x15x15" (33x38x38cm)
<b>N</b> 33-02005YF	20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm)	100	14x9x19" (40x23x48cm)
<b>S</b> 33-02005YF-S	20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm)	80	13x15x15" (33x38x38cm)
<b>N</b> 33-02007	20-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-02007-S	20-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-02007Y	20-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-02007Y-S	20-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-02010	20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x5" (5x13cm)	100	14x9x19" (40x23x48cm)
<b>S</b> 33-02010-S	20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x5" (5x13cm)	80	13x15x15" (33x38x38cm)
<b>N</b> 33-02010F	20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x5" (5x13cm)	100	13x15x15" (33x38x38cm)
<b>S</b> 33-02010F-S	20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x1x5" (5x3x13cm)	80	13x15x15" (33x38x38cm)
<b>N</b> 33-02010Y	20-Count, Foam Block, Adhesive, Blade Remover, Yellow, 2x5" (5x13cm)	100	15x9x20" (36x22x50cm)
<b>S</b> 33-02010Y-S	20-Count, Foam Block, Adhesive, Blade Remover, Yellow, 2x5" (5x13cm)	80	13x15x15" (33x38x38cm)
<b>N</b> 33-02080	20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-02080-S	20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-02080Y	20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-02080Y-S	20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-02081	20-Count, Foam Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-02081-S	20-Count, Foam Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-02081Y	20-Count, Foam Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-02081Y-S	20-Count, Foam Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-03007	30-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-03007-S	30-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-03007Y	30-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-03007Y-S	30-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-03080	30-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 4x5" (10x13cm)	100	17x6x21" (43x15x53cm)
<b>S</b> 33-03080-S	30-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 4x5" (10x13cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 33-03080Y	30-Count, Foam Strip, Magnet Strip, Adhesive, Blade Remover, Yellow, 5x4" (12x9cm)	100	21x17x6" (54x44x14cm)
<b>N</b> 33-03085	30-Count, Foam Strip, Magnet, Blade Remover, Red, 4x5" (10x13cm)	100	17x6x21" (43x15x53cm)
<b>S</b> 33-03085-S	30-Count, Foam Strip, Magnet, Blade Remover, Red, 4x5" (10x13cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 33-03085Y	30-Count, Foam Strip, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	100	17x6x21" (43x15x53cm)
<b>S</b> 33-03085Y-S	30-Count, Foam Strip, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 33-03088	30-Count, Foam Strips, Blade Remover, Red, 2x3" (5x7cm)	100	12x6x18" (30x15x46cm)

# Surgical Sharps Safety

## SECURE Needle Counters

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
<b>S</b> 33-03088-S	30-Count, Foam Strips, Blade Remover, Red, 2x3" (5x7cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-03088Y	30-Count, Foam Strips, Blade Remover, Yellow, 4x5" (10x13cm)	100	12x6x18" (30x15x46cm)
<b>S</b> 33-03088Y-S	30-Count, Foam Strips, Blade Remover, Yellow, 4x5" (10x13cm)	100	13x15x15" (33x38x38cm)
<b>N</b> 33-03505Y	Counter, Foam Block, Magnet, Blade Remover 4x2" (11x6cm)	100	20x14x9" (50x36x22cm)
<b>N</b> 33-04010	40-Count, Foam Block, Adhesive, Blade Remover, Red, 4x5" (10x13cm)	100	16x15x18" (41x38x46cm)
<b>S</b> 33-04010-S	40-Count, Foam Block, Adhesive, Blade Remover, Red, 4x5" (10x13cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 33-04010Y	40-Count, Foam Block, Adhesive, Blade Remover, Yellow, 4x5" (10x13cm)	100	16x15x18" (41x38x46cm)
<b>S</b> 33-04010Y-S	40-Count, Foam Block, Adhesive, Blade Remover, Yellow, 4x5" (10x13cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 33-04037	40-Count, Foam Block, Magnet, Blade Remover, Red, 4x5" (10x13cm)	100	16x15x18" (41x38x46cm)
<b>S</b> 33-04037-S	40-Count, Foam Block, Magnet, Blade Remover, Red, 4x5" (10x13cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 33-04037Y	40-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	100	16x15x18" (41x38x46cm)
<b>S</b> 33-04037Y-S	40-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 33-04037YF	40-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	100	16x15x18" (41x38x46cm)
<b>N</b> 33-06007	60-Count, Double Magnet, Blade Remover, Red, 4x5" (10x13cm)	100	17x6x21" (43x15x53cm)
<b>S</b> 33-06007-S	60-Count, Double Magnet, Blade Remover, Red, 4x5" (10x13cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 33-06007Y	60-Count, Double Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	100	17x6x21" (43x15x53cm)
<b>S</b> 33-06007Y-S	60-Count, Double Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 33-07037	10-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	100	16x15x18" (40x38x46cm)
<b>S</b> 33-07037-S	10-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 33-07037Y	10-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm)	100	16x15x18" (40x38x46cm)
<b>S</b> 33-07037Y-S	70-Count, Foam Block, Magnet w/count, Blade Remover, Yellow, 4x5" (10x13cm)	50	15x11x25" (38x28x64cm)

### Materials

- Colored Boxes: HIPS Plastic
- Clear Lids: Acrylonitrile Butadiene Styrene (ABS)
- Foam: Polyethylene (PE)
- Magnet: Iron and Polyethylene

### Product Attributes

- Specially selected foam does not shrink during sterilization
- Slide locks and clear lids available. Optional blade removers

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
 www.alleset.com  
 SM-TDS-003 2024  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
 www.alleset.com  
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
 Westervoortsedijk 60,  
 6827 AT Arnhem,  
 The Netherlands  
 www.emergobyul.com

