

Surgical Sharps Safety

SECURE Needle Counters

| CODE | DESCRIPTION | CASE QTY | CASE DIMENSIONS |
|-----------------------|--|----------|------------------------|
| N 33-01037 | 10-Count, Foam Block, Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 15x13x7" (38x33x18cm) |
| S 33-01037-S | 10-Count, Foam Block, Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 15x13x7" (38x33x18cm) |
| N 33-01037Y | 10-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 15x13x7" (38x33x18cm) |
| S 33-01037Y-S | 10-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 15x13x7" (38x33x18cm) |
| N 33-01580 | 15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x5" (5x13cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-01580-S | 15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x5" (5x13cm) | 100 | 13x10x15" (33x25x38cm) |
| N 33-01580Y | 15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x5" (5x13cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-01580Y-S | 15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow, 2x5" (5x13cm) | 100 | 13x10x15" (33x25x38cm) |
| N 33-01585 | 15-Count, Foam Strips, Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-01585-S | 15-Count, Foam Strips, Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 13x10x15" (33x25x38cm) |
| N 33-01585Y | 15-Count, Foam Strips, Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-01585Y-S | 15-Count, Foam Strips, Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 13x10x15" (33x25x38cm) |
| N 33-02005 | 20-Count, Foam Block, Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 14x9x20" (36x22x50cm) |
| S 33-02005-S | 20-Count, Foam Block, Magnet, Blade Remover, Red, 2x5" (5x13cm) | 80 | 13x15x15" (33x38x38cm) |
| N 33-02005Y | 20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm) | 100 | 14x9x20" (36x22x50cm) |
| S 33-02005Y-S | 20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm) | 80 | 13x15x15" (33x38x38cm) |
| N 33-02005YF | 20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm) | 100 | 14x9x19" (40x23x48cm) |
| S 33-02005YF-S | 20-Count, Foam Block, Magnet, Blade Remover, Yellow, 2x5" (5x13cm) | 80 | 13x15x15" (33x38x38cm) |
| N 33-02007 | 20-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-02007-S | 20-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-02007Y | 20-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-02007Y-S | 20-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-02010 | 20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x5" (5x13cm) | 100 | 14x9x19" (40x23x48cm) |
| S 33-02010-S | 20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x5" (5x13cm) | 80 | 13x15x15" (33x38x38cm) |
| N 33-02010F | 20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x5" (5x13cm) | 100 | 13x15x15" (33x38x38cm) |
| S 33-02010F-S | 20-Count, Foam Block, Adhesive, Blade Remover, Red, 2x1x5" (5x3x13cm) | 80 | 13x15x15" (33x38x38cm) |
| N 33-02010Y | 20-Count, Foam Block, Adhesive, Blade Remover, Yellow, 2x5" (5x13cm) | 100 | 15x9x20" (36x22x50cm) |
| S 33-02010Y-S | 20-Count, Foam Block, Adhesive, Blade Remover, Yellow, 2x5" (5x13cm) | 80 | 13x15x15" (33x38x38cm) |
| N 33-02080 | 20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-02080-S | 20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-02080Y | 20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-02080Y-S | 20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-02081 | 20-Count, Foam Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-02081-S | 20-Count, Foam Strip, Adhesive, Blade Remover, Red, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-02081Y | 20-Count, Foam Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-02081Y-S | 20-Count, Foam Strip, Adhesive, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-03007 | 30-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-03007-S | 30-Count, Double Magnet, Blade Remover, Red, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-03007Y | 30-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-03007Y-S | 30-Count, Double Magnet, Blade Remover, Yellow, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-03080 | 30-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 4x5" (10x13cm) | 100 | 17x6x21" (43x15x53cm) |
| S 33-03080-S | 30-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Red, 4x5" (10x13cm) | 108 | 15x11x25" (38x28x64cm) |
| N 33-03080Y | 30-Count, Foam Strip, Magnet Strip, Adhesive, Blade Remover, Yellow, 5x4" (12x9cm) | 100 | 21x17x6" (54x44x14cm) |
| N 33-03085 | 30-Count, Foam Strip, Magnet, Blade Remover, Red, 4x5" (10x13cm) | 100 | 17x6x21" (43x15x53cm) |
| S 33-03085-S | 30-Count, Foam Strip, Magnet, Blade Remover, Red, 4x5" (10x13cm) | 108 | 15x11x25" (38x28x64cm) |
| N 33-03085Y | 30-Count, Foam Strip, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 17x6x21" (43x15x53cm) |
| S 33-03085Y-S | 30-Count, Foam Strip, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 108 | 15x11x25" (38x28x64cm) |
| N 33-03088 | 30-Count, Foam Strips, Blade Remover, Red, 2x3" (5x7cm) | 100 | 12x6x18" (30x15x46cm) |

Surgical Sharps Safety

SECURE Needle Counters

| CODE | DESCRIPTION | CASE QTY | CASE DIMENSIONS |
|----------------------|---|----------|------------------------|
| S 33-03088-S | 30-Count, Foam Strips, Blade Remover, Red, 2x3" (5x7cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-03088Y | 30-Count, Foam Strips, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 12x6x18" (30x15x46cm) |
| S 33-03088Y-S | 30-Count, Foam Strips, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 13x15x15" (33x38x38cm) |
| N 33-03505Y | Counter, Foam Block, Magnet, Blade Remover 4x2" (11x6cm) | 100 | 20x14x9" (50x36x22cm) |
| N 33-04010 | 40-Count, Foam Block, Adhesive, Blade Remover, Red, 4x5" (10x13cm) | 100 | 16x15x18" (41x38x46cm) |
| S 33-04010-S | 40-Count, Foam Block, Adhesive, Blade Remover, Red, 4x5" (10x13cm) | 50 | 15x11x25" (38x28x64cm) |
| N 33-04010Y | 40-Count, Foam Block, Adhesive, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 16x15x18" (41x38x46cm) |
| S 33-04010Y-S | 40-Count, Foam Block, Adhesive, Blade Remover, Yellow, 4x5" (10x13cm) | 50 | 15x11x25" (38x28x64cm) |
| N 33-04037 | 40-Count, Foam Block, Magnet, Blade Remover, Red, 4x5" (10x13cm) | 100 | 16x15x18" (41x38x46cm) |
| S 33-04037-S | 40-Count, Foam Block, Magnet, Blade Remover, Red, 4x5" (10x13cm) | 50 | 15x11x25" (38x28x64cm) |
| N 33-04037Y | 40-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 16x15x18" (41x38x46cm) |
| S 33-04037Y-S | 40-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 50 | 15x11x25" (38x28x64cm) |
| N 33-04037YF | 40-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 16x15x18" (41x38x46cm) |
| N 33-06007 | 60-Count, Double Magnet, Blade Remover, Red, 4x5" (10x13cm) | 100 | 17x6x21" (43x15x53cm) |
| S 33-06007-S | 60-Count, Double Magnet, Blade Remover, Red, 4x5" (10x13cm) | 108 | 15x11x25" (38x28x64cm) |
| N 33-06007Y | 60-Count, Double Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 17x6x21" (43x15x53cm) |
| S 33-06007Y-S | 60-Count, Double Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 108 | 15x11x25" (38x28x64cm) |
| N 33-07037 | 10-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 16x15x18" (40x38x46cm) |
| S 33-07037-S | 10-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 50 | 15x11x25" (38x28x64cm) |
| N 33-07037Y | 10-Count, Foam Block, Magnet, Blade Remover, Yellow, 4x5" (10x13cm) | 100 | 16x15x18" (40x38x46cm) |
| S 33-07037Y-S | 70-Count, Foam Block, Magnet w/count, Blade Remover, Yellow, 4x5" (10x13cm) | 50 | 15x11x25" (38x28x64cm) |

Materials

- Colored Boxes: HIPS Plastic
- Clear Lids: Acrylonitrile Butadiene Styrene (ABS)
- Foam: Polyethylene (PE)
- Magnet: Iron and Polyethylene

Product Attributes

Specially selected foam does not shrink during sterilization. Slide locks and clear lids available. Optional blade removers.

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Establishes at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.



Product Certifications


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
 www.alleset.com
 SM-TDS-002 2022
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
 www.alleset.com
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 www.emergobyul.com

