

















































# Surgical Sharps Safety

## SELECT Needle Counters

CODE	DESCRIPTION (AxB)	CASE QTY	CASE DIMENSIONS
 30-01001	10-Count, Magnet, Clear Top, Red Base, 2x4" (5.5x11cm)	100	10x7x13" (26x18x33cm)
 30-01001-S	10-Count, Magnet, Clear Top, Red Base, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01001Y	10-Count, Magnet, Clear Top, Yellow Base, 2x4" (5.5x11cm)	100	10x7x13" (26x18x33cm)
 30-01001Y-S	10-Count, Magnet, Clear Top, Yellow Base, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01010	10-Count, Foam Block, Red, 2x3" (5.5x7.5cm)	100	13x7x15" (33x18x38cm)
 30-01010-S	10-Count, Foam Block, Red, 2x3" (5.5x7.5cm)	100	13x13x15" (33x33x38cm)
 30-01010Y	10-Count, Foam Block, Yellow, 2x3" (5.5x7.5cm)	100	13x7x15" (33x18x38cm)
 30-01010Y-S	10-Count, Foam Block, Yellow, 2x3" (5.5x7.5cm)	100	13x13x15" (33x33x38cm)
 30-01037	10-Count, Foam Block, Magnet, Red, 2x3" (5.5x7.5cm)	100	13x7x15" (33x18x38cm)
 30-01037-S	10-Count, Foam Block, Magnet, Red, 2x3" (5.5x7.5cm)	100	13x13x15" (33x33x38cm)
 30-01037Y	10-Count, Foam Block, Magnet, Yellow, 2x3" (5.5x7.5cm)	100	13x7x15" (33x18x38cm)
 30-01037Y-S	10-Count, Foam Block, Magnet, Yellow, 2x3" (5.5x7.5cm)	100	13x13x15" (33x33x38cm)
 30-01037YF	10-Count, Foam Block, Magnet, Yellow, 2x3" (5.5x7.5cm)	100	13x7x15" (33x18x38cm)
 30-01037YF-S	10-Count, Foam Block, Magnet, Yellow, 2x3" (5.5x7.5cm)	100	13x13x15" (33x33x38cm)
 30-01514	15-Count, Magnet, Clear Top, Red Base, 2x4" (5.5x11cm)	168	14x6x18" (36x15x46cm)
 30-01514-S	15-Count, Magnet, Clear Top, Red Base, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01514Y	15-Count, Magnet, Clear Top, Yellow Base, 2x4" (5.5x11cm)	168	14x6x18" (36x15x46cm)
 30-01514Y-S	15-Count, Magnet, Clear Top, Yellow Base, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01518	15-Count, Foam Strips, Clear Top, Red Base, 2x4" (5.5x11cm)	96	9x5x18" (24x12x46cm)
 30-01518-S	15-Count, Foam Strips, Clear Top, Red Base, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01518Y	15-Count, Foam Strips, Clear Top, Yellow Base, 2x4" (5.5x11cm)	96	9x5x18" (24x12x46cm)
 30-01518Y-S	15-Count, Foam Strips, Clear Top, Yellow Base, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01585	15-Count, Foam Strips, Magnet, Red, 2x4" (5.5x11cm)	168	14x6x18" (36x14x46cm)
 30-01585-S	15-Count, Foam Strips, Magnet, Red, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-01585Y	15-Count, Foam Strips, Magnet, Yellow, 2x4" (5.5x11cm)	168	14x6x18" (36x14x46cm)
 30-01585Y-S	15-Count, Foam Strips, Magnet, Yellow, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-02005	20-Count, Foam Block, Magnet, Red, 2x4" (5.5x11cm)	100	10x13x14" (26x33x36cm)
 30-02005-S	20-Count, Foam Block, Magnet, Red, 2x4" (5.5x11cm)	100	13x13x15" (33x33x38cm)
 30-02005Y	20-Count, Foam Block, Magnet, Yellow, 2x4" (5.5x11cm)	100	10x13x14" (26x33x36cm)
 30-02005Y-S	20-Count, Foam Block, Magnet, Yellow, 2x4" (5.5x11cm)	100	13x13x15" (33x33x38cm)
 30-02010	20-Count, Foam Block, Adhesive, Red, 2x4" (5.5x11cm)	100	10x13x14" (26x33x36cm)
 30-02010-S	20-Count, Foam Block, Adhesive, Red, 2x4" (5.5x11cm)	100	13x13x15" (33x33x38cm)
 30-02010Y	20-Count, Foam Block, Adhesive, Yellow, 2x4" (5.5x11cm)	100	10x13x14" (26x33x36cm)
 30-02010Y-S	20-Count, Foam Block, Adhesive, Yellow, 2x4" (5.5x11cm)	100	13x13x15" (33x33x38cm)
 30-02055	20-Count, Double Magnet, Red, 2x4" (5.5x11cm)	168	14x6x18" (36x15x46cm)
 30-02055-S	20-Count, Double Magnet, Red, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-02055Y	20-Count, Double Magnet, Yellow, 2x4" (5.5x11cm)	168	14x6x18" (36x15x46cm)
 30-02055Y-S	20-Count, Double Magnet, Yellow, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-03007	30-Count, Double Magnet, Red, 2x4" (5.5x11cm)	100	10x7x13" (26x18x33cm)
 30-03007-S	30-Count, Double Magnet, Red, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-03007Y	30-Count, Double Magnet, Yellow, 2x4" (5.5x11cm)	100	10x7x13" (26x18x33cm)
 30-03007Y-S	30-Count, Double Magnet, Yellow, 2x4" (5.5x11cm)	100	13x10x15" (33x25x38cm)
 30-03015	30-Count, Magnet, Clear Top, Red Base, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
 30-03015-S	30-Count, Magnet, Clear Top, Red Base, 4x5" (10x12cm)	108	15x10x24" (37x26x62cm)
 30-03015Y	30-Count, Magnet, Clear Top, Yellow Base, 4x5" (10x12cm)	64	12x6x18" (30x15x46cm)
 30-03015Y-S	30-Count, Magnet, Clear Top, Yellow Base, 4x5" (10x12cm)	108	15x10x24" (37x26x62cm)
 30-03085	30-Count, Foam Strips, Magnet, Red, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
 30-03085-S	30-Count, Foam Strips, Magnet, Red, 4x5" (10x12cm)	108	15x10x24" (37x26x62cm)

# Surgical Sharps Safety

## SELECT Needle Counters

CODE	DESCRIPTION (AxB)	CASE QTY	CASE DIMENSIONS
<b>N</b> 30-03085Y	30-Count, Foam Strips, Magnet, Yellow, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
<b>S</b> 30-03085Y-S	30-Count, Foam Strips, Magnet, Yellow, 4x5" (10x12cm)	108	15x10x24" (37x26x62cm)
<b>N</b> 30-03088	30-Count, Foam Strips, Red, 2x4" (5.5x7.5cm)	96	10x7x13" (26x18x33cm)
<b>S</b> 30-03088-S	30-Count, Foam Strips, Red, 2x4" (5.5x7.5cm)	100	13x10x15" (33x25x38cm)
<b>N</b> 30-03088Y	30-Count, Foam Strips, Yellow, 2x4" (5.5x7.5cm)	96	10x7x13" (26x18x33cm)
<b>S</b> 30-03088Y-S	30-Count, Foam Strips, Yellow, 2x4" (5.5x7.5cm)	100	13x10x15" (33x25x38cm)
<b>N</b> 30-04027	40-Count, Foam Block, Adhesive, Red, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-04027-S	40-Count, Foam Block, Adhesive, Red, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-04027Y	40-Count, Foam Block, Adhesive, Yellow, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-04027Y-S	40-Count, Foam Block, Adhesive, Yellow, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-04037	40-Count, Foam Block, Magnet, Red, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-04037-S	40-Count, Foam Block, Magnet, Red, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-04037Y	40-Count, Foam Block, Magnet, Yellow, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-04037Y-S	40-Count, Foam Block, Magnet, Yellow, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-06055	60-Count, Double Magnet, Red, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
<b>S</b> 30-06055-S	60-Count, Double Magnet, Red, 4x5" (10x12cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 30-06055Y	60-Count, Double Magnet, Yellow, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
<b>S</b> 30-06055Y-S	60-Count, Double Magnet, Yellow, 4x5" (10x12cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 30-06088	60-Count, Foam Strips, Red, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
<b>S</b> 30-06088-S	60-Count, Foam Strips, Red, 4x5" (10x12cm)	108	15x11x25" (38x28x64cm)
<b>N</b> 30-06088Y	60-Count, Foam Strips, Yellow, 4x5" (10x12cm)	64	9x6x17" (24x16x44cm)
<b>S</b> 30-06088Y-S	60-Count, Foam Strips, Yellow, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-08077	80-Count, Double Foam Block, Red, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-08077-S	80-Count, Double Foam Block, Red, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-08077Y	80-Count, Double Foam Block, Yellow, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-08077Y-S	80-Count, Double Foam Block, Yellow, 4x5" (10x12cm)	50	15x11x25" (38x28x64cm)
<b>N</b> 30-10055	100-Count, Double Magnet, Foam Block, Red, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-10055-S	100-Count, Double Magnet, Foam Block, Red, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>N</b> 30-10055Y	100-Count, Double Magnet, Foam Block, Yellow, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)
<b>S</b> 30-10055Y-S	100-Count, Double Magnet, Foam Block, Yellow, 4x5" (10x12cm)	64	16x6x20" (40x16x50cm)

# Surgical Sharps Safety

## SELECT Needle Counters

### Materials

- Colored Boxes: HIPS Plastic
- Clear Lids: Acrylonitrile Butadiene Styrene (ABS)
- Foam: Polyethylene (PE)
- Magnet: Iron and Polyethylene

### Product Attributes

- Specially selected foam does not shrink during sterilization
- Slide locks and clear lids available
- Optional blade removers

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
[www.alleset.com](http://www.alleset.com)  
 SM-TDS-003 2023  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
 Westervoortsedijk 60,  
 6827 AT Arnhem,  
 The Netherlands  
[www.alleset.com](http://www.alleset.com)