



























# Orthopedic Dressing

## AlleFlex™ Cohesive Bandages

| CODE  | DESCRIPTION (WxL)   | CASE QTY | CASE DIMENSIONS      |
|---|---|----------|----------------------|
|  86-10CA   | AlleFlex™ Cohesive Bandages, Camo, 1x180" (2.5x457cm)         | 48       | 10x7x6" (26x18x14cm) |
|  86-10PR   | AlleFlex™ Cohesive Bandages, Prism, 1x180" (2.5x457cm)        | 48       | 10x7x6" (26x18x14cm) |
|  86-10PW   | AlleFlex™ Cohesive Bandages, Paw Print, 1x180" (2.5x457cm)    | 48       | 10x7x6" (26x18x14cm) |
|  86-10TN   | AlleFlex™ Cohesive Bandages, Tan, 1x180" (2.5x457cm)          | 48       | 10x7x6" (26x18x14cm) |
|  86-15CA   | AlleFlex™ Cohesive Bandages, Camo, 1.5x180" (3.8x457cm)       | 36       | 12x7x6" (30x18x14cm) |
|  86-15CK   | AlleFlex™ Cohesive Bandages, Checkered, 1.5x180" (3.8 x457cm) | 36       | 12x7x6" (30x18x14cm) |
|  86-15PR   | AlleFlex™ Cohesive Bandages, Prism, 1.5x180" (3.8 x457cm)     | 36       | 12x7x6" (30x18x14cm) |
|  86-15PW   | AlleFlex™ Cohesive Bandages, Paw Print, 1.5x180" (3.8x457cm)  | 36       | 12x7x6" (30x18x14cm) |
|  86-15TN   | AlleFlex™ Cohesive Bandages, Tan, 1.5x180" (3.8x457cm)        | 36       | 12x7x6" (30x18x14cm) |
|  86-20CA   | AlleFlex™ Cohesive Bandages, Camo, 2x180" (5x457cm)           | 36       | 15x7x6" (37x18x14cm) |
|  86-20CK   | AlleFlex™ Cohesive Bandages, Checkered, 2x180" (5x457cm)      | 36       | 15x7x6" (37x18x14cm) |
|  86-20DB   | AlleFlex™ Cohesive Bandages, Dark Blue, 2x180" (5x457cm)      | 36       | 15x7x6" (37x18x14cm) |
|  86-20KP   | AlleFlex™ Cohesive Bandages, Kids, 2x180" (5x457cm)           | 36       | 15x7x6" (37x18x14cm) |
|  86-20PR   | AlleFlex™ Cohesive Bandages, Prism, 2x180" (5x457cm)          | 36       | 15x7x6" (37x18x14cm) |
|  86-20PW   | AlleFlex™ Cohesive Bandages, Paw Print, 2x180" (5x457cm)      | 36       | 15x7x6" (37x18x14cm) |
|  86-20TN   | AlleFlex™ Cohesive Bandages, Tan, 2x180" (5x457cm)            | 36       | 15x7x6" (37x18x14cm) |
|  86-30CA   | AlleFlex™ Cohesive Bandages, Camo, 3x180" (7.6x457cm)         | 24       | 15x7x6" (37x18x14cm) |
|  86-30CK  | AlleFlex™ Cohesive Bandages, Checkered, 3x180" (7.6x457cm)    | 24       | 15x7x6" (37x18x14cm) |
|  86-30PR | AlleFlex™ Cohesive Bandages, Prism, 3x180" (7.6x457cm)        | 24       | 15x7x6" (37x18x14cm) |
|  86-30PW | AlleFlex™ Cohesive Bandages, Paw Print, 3x180" (7.6x457cm)    | 24       | 15x7x6" (37x18x14cm) |
|  86-30TN | AlleFlex™ Cohesive Bandages, Tan, 3x180" (7.6x457cm)          | 24       | 15x7x6" (37x18x14cm) |
|  86-40CA | AlleFlex™ Cohesive Bandages, Camo, 4x180" (10x457cm)          | 18       | 9x7x8" (24x18x20cm)  |
|  86-40PR | AlleFlex™ Cohesive Bandages, Prism, 4x180" (10x457cm)         | 18       | 9x7x8" (24x18x20cm)  |
|  86-40PW | AlleFlex™ Cohesive Bandages, Paw Print, 4x180" (10x457cm)     | 18       | 9x7x8" (24x18x20cm)  |
|  86-40TN | AlleFlex™ Cohesive Bandages, Tan, 4x180" (10x457cm)           | 18       | 9x7x8" (24x18x20cm)  |
|  86-60TN | AlleFlex™ Cohesive Bandages, Tan, 6x180" (15x457cm)           | 18       | 13x7x8" (34x18x20cm) |

### Materials

- Non-Woven Material and Elastic Fibers with Proprietary Self-Adherent Coating

### Product Attributes


- Softspun fabric for easy, quick tearing
- Industry leading self-adherence

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.



# Orthopedic Dressing

## AlleFlex™ Cohesive Bandages

### Product Certifications/Standards

---

- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/ Registrations

---

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

---

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
[www.alleset.com](http://www.alleset.com)  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
[www.emergobyul.com](http://www.emergobyul.com)

