

Neonatal and Infant Baby Blankets & Caps

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
25-1336	Baby Cap, 6x4" (15x10cm)	1400	24x16x16" (62x40x41cm)
25-1340	Baby Blanket, Pink/Blue Stripes, 30x41" (75x103cm)	100	9x17x21" (24x42x52cm)
25-1341	Baby Blanket, Pink/Blue Babies, 30x41" (75x103cm)	100	9x17x21" (24x42x52cm)

Materials

- 25-1340, 25-1341: Cotton
- 25-1336: Polyester

Product Attributes

- Premium soft-napped fabric for enhanced softness and comfort
- Disposable materials with same feel as standard blankets

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions



25-1336



25-1340



25-1341



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003-2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

