

# Surgical Drapes

## ECO<sup>2</sup> Back Table Covers

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
E80-13129	ECO <sup>2</sup> Zone Reinforced Back Table Cover, 44x44" (112x112cm)	60	22x12x12" (56x31x31cm)
E80-13129-HB	ECO <sup>2</sup> Zone Reinforced Back Table Cover, Handi-Bin, 44x44" (112x112cm)	640	40x24x18" (102x61x45cm)

### Materials

- ECO<sup>2</sup> Hydrophilic Spunbond Polypropylene (PP)
- ECO<sup>2</sup> Blue Polyethylene (PE) Film

### Product Attributes

- Certain components of this product incorporate ECO<sup>2</sup> sustainable materials which enhance product biodegradability
- Absorbent reinforcement material

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

