

Surgical Drapes

ECO² Back Table Covers

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
E80-13132	ECO ² Fully Reinforced Back Table cover, 59x94" (150x240cm)	32	18x15x11" (45x37x28cm)
E80-13132-S	ECO ² Fully Reinforced Back Table cover, 59x94" (150x240cm)	26	23x15x16" (59x39x40cm)
E80-13132-HB	ECO ² Fully Reinforced Back Table cover, Handi-Bin, 59x94" (150x240cm)	320	46x30x31" (116x75x78cm)

Materials

- ECO² ONEZONE[®] Bi-Laminate Nonwoven Material (PP/PE)

Product Attributes

- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability
- Absorbent top layer with impermeable film bottom layer for maximum protection

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 3 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

