

Surgical Drapes

Back Table Covers

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
80-13166-HB	Back Table Cover, Handi-Bin, 79x90" (200x228cm)	192	47x31x16" (119x79x40cm)

Materials

- Hydrophilic Spunbond Polypropylene (PP)
- Blue Polyethylene (PE) Film

Product Attributes

- Zone reinforced

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleaset.

Shelf Life

Non-Sterile: Establishes at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions




Alleaset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleaset.com

SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleaset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com