

Surgical Drapes

Back Table Covers

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
80-13322	Back Table Cover, 60x90" (152x229cm)	14	18x13x11" (45x34x29cm)
80-13322-S	Back Table Cover, 60x90" (152x229cm)	26	20x14x20" (52x35x52cm)
80-13322-HB	Back Table Cover, Handi-Bin, 60x90" (152x229cm)	140	39x31x17" (100x80x44cm)

Materials

- Blue Polyethylene (PE) Film
- HydraGuard® Absorbent Reinforcement

Product Attributes

- Heavy-duty film
- Absorbent reinforcement material for fluid management

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleset.com
 SM-TDS-003 2023
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.
 Westervoortsedijk 60, 6827
 AT Arnhem,
 The Netherlands.
www.emergobyul.com

