

Surgical Accessories

Bulb Syringes

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
25-6700	Syringe Bulb, 60ml (2oz)	50	12x12x12" (30x30x30cm)
25-6700-S	Syringe Bulb, 60ml (2oz)	TBD	TBD
25-6702	Syringe Bulb, 60ml (2oz), Blue	50	12x12x12" (30x30x30cm)
25-6702-S	Syringe Bulb, 60ml (2oz), Blue	200	18x13x15" (46x33x37cm)
25-6703	Syringe Bulb, 89ml (3oz), Green	50	12x12x12" (30x30x30cm)
25-6703-S	Syringe Bulb, 89ml (3oz), Green	150	19x16x10" (49x41x26cm)

Materials

- Polyvinyl Chloride (PVC)
- Polypropylene (PP)

Product Attributes

- Pliable bulb for responsive, controllable suction
- Finger flanges to help prevent roll

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.



Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards

- CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 5

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2024
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands.
www.emergobyul.com

