

Equipment Covers

Endo Camera

CODE	DESCRIPTION (WxL)	CASE QTY	CASE DIMENSIONS
N 44-0596EV	Camera Cover, Elastic Tip, w/ Cardboard Sleeve Insert, 5x96" (13x244cm)	100	18x12x6" (46x28x14cm)
S 44-0596EV-S	Camera Cover, Elastic Tip, w/ Cardboard Sleeve Insert, 5x96" (13x244cm)	50	18x12x6" (46x28x14cm)
N 44-0796EV	Camera Cover, Elastic Tip, w/ Cardboard Sleeve Insert, 7x96" (18x244cm)	100	11x9x10" (28x23x25cm)
S 44-0796EV-S	Camera Cover, Elastic Tip, w/ Cardboard Sleeve Insert, 7x96" (18x244cm)	50	11x9x10" (28x23x25cm)
N 44-0996EV	Camera Cover, Elastic Tip, w/ Cardboard Sleeve Insert, 9x96" (23x244cm)	100	15x12x12" (38x30.5x30cm)
S 44-0996EV-S	Camera Cover, Elastic Tip, w/ Cardboard Sleeve Insert, 9x96" (23x244cm)	50	15x12x12" (38x30.5x30cm)

Materials

- Clear Polyethylene (PE) Film
- Conformable Fenestration Material
- Cardboard

Product Attributes

- Perforated conformable film tip
- Adhesive polyethylene securement strap
- Cardboard sleeve opening

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleaset.com
 SM-TDS-003 2024
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleaset.com
 US FDA Establishment Registration #3004911384

EC REP Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.emergobyul.com