

Temperature Maintenance

Cozee™ Caregiver Attire

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
68-480311	Pediatric Jacket, Green Sleeve, Silver Body	XS	50	13x9x13" (34x24x32cm)
68-480321	Pediatric Jacket, Green Sleeve, Silver Body	S	50	13x9x13" (34x24x32cm)
68-480331	Pediatric Jacket, Green Sleeve, Silver Body	M	50	18x12x8" (46x30x21cm)
68-480341	Pediatric Jacket, Green Sleeve, Silver Body	L	50	18x12x8" (46x30x21cm)
68-480362	Jacket, Green w/ Silver Arms	L	25	17x13x10" (42x32x26cm)
68-480352	Jacket, Green w/ Silver Arms	M	25	17x15x11" (42x32x26cm)
68-480351	Jacket, All Silver Green Inside	M	25	17x13x10" (42x32x26cm)
68-480361	Jacket, All Silver Green Inside	L	25	17x15x11" (42x37x28cm)
68-480371	Jacket, All Silver Green Inside	XL	25	17x13x10" (42x32x26cm)
68-480372	Jacket, Green w/ Silver Arms	XL	25	17x13x10" (42x32x26cm)
68-480841	Vest, Silver, Green Inside	S	25	13x11x8" (32x28x21.5cm)
68-480851	Vest, Silver, Green Inside	M	25	12x10x9" (31x25x22cm)
68-480861	Vest, Silver, Green Inside	L	25	13x11x8" (32x28x21.5cm)
68-480871	Vest, Silver, Green Inside	XL	25	13x11x8" (32x28x21.5cm)
68-480711	Pediatric Pant	XS	50	18x9x6" (46x24x16cm)
68-480721	Pediatric Pant	S	50	18x9x6" (46x24x16cm)
68-480731	Pediatric Pant	M	50	20x12x7" (51x30x18cm)
68-480741	Pediatric Pant	L	50	20x12x11" (52x30x28cm)

Materials

- Non-Woven Green Polypropylene (PP)
- Aluminum Laminate on Polyethylene (PE) film

Product Attributes

- High thermal retention
- Non-conductive

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Not Applicable

Product Certifications


- CE certified as class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/ Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions



68-480311, 68-480321,
68-480331 68-480341



68-480841 68-480851
68-480861 68-480871



68-480711, 68-480721,
68-480731, 68-480741



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

