

Temperature Maintenance

Cozee™ Patient Attire

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
68-480221	Cap, Neonatal w/ Stockinette Band	-	100	13x7x5" (33x18x12cm)
68-480231	Cap, Infant	-	100	12x9x10" (31x23x25cm)
68-480241	Cap, Child	-	50	13x7x5" (33x18x12cm)
68-480251	Cap, Youth	-	50	12x11x10" (31x28x26cm)
68-480261	Cap, Adult	-	100	25x14x9" (63x35x23cm)
68-480441	Pulse Oximeter Mitt, Infant	S	100	15x7x6" (37x18x14cm)
68-480451	Pulse Oximeter Mitt, Infant	M	100	12x11x8" (31x27x20cm)
68-480461	Pulse Oximeter Mitt, Infant	L	50	13x7x5" (33x18x12cm)
68-480562	Patient Gown, Silver Front w/ Green Arms	-	50	21x15x14" (54x39x35.5cm)
68-480651	Foot Warmers	-	50	13x10x6" (32x26x15cm)
68-480751	Adjustable Leg warmers	-	100	19x17x11" (48x43x27cm)

Materials

- Non-Woven Green Polypropylene (PP)
- Aluminum Laminate on Polyethylene (PE) film

Product Attributes

- High thermal retention
- Non-conductive

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Not Applicable

Product Certifications


- CE certified as class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

Storage


- General warehouse conditions



68-480221, 68-480231

68-480441, 68-480451
68-48461

68-480651



68-480562

68-480751



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

