

Temperature Maintenance

Cozee™ Warming Blankets and Wraps

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
68-4115137	Warming Blanket, Non-Reflective, 51x37" (130x95cm)	S	40	24x15x14" (60x37x35cm)
68-4115175	Warming Blanket, Non-Reflective, 51x75" (130x190cm)	L	20	24x15x14" (60x37x35cm)
68-480141	Warming Blanket, Reflective, 35x40" (88x100cm)	S	60	12x9x10" (31x23x25cm)
68-480151	Warming Blanket, Reflective, 48x48" (120x120cm)	M	50	18x12x7" (46x30x18cm)
68-480161	Warming Blanket, Reflective, 48x84" (122x213cm)	L	25	14x9x13" (34x24x32cm)
68-480101	Papoose Infant Wrap, Reflective	-	50	12x10x9" (31x25x22cm)
68-480901	Rescue Wrap, Reflective, Full Enclosure Bag	-	10	17x15x11" (42x37x28cm)

Materials

Warming Blankets, Non-Reflective

- Non-Woven Polypropylene (PP)
- Polyethylene Terephthalate (PET)

Warming Blankets and Wraps, Reflective

- Non-Woven Green Polypropylene (PP)
- Aluminum Laminate on Polyethylene (PE) Film
- Nylon Hook and Loop

Product Attributes

Non-reflective	Reflective	Sharing traits
<ul style="list-style-type: none"> • Three-layer warmth • 20% better insulation than cotton blankets 	<ul style="list-style-type: none"> • Reflective • Non-conductive • Non-linting 	<ul style="list-style-type: none"> • High thermal retention • Individually shrinked wrapped

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Not Applicable

Product Certifications


• CE certified as class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

• EN ISO 13485:2016

Storage

• General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com