

Temperature Maintenance

Cozee™ Kit

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
68-480951	Kit, Blanket w/ Cap and 2 Leg Warmers	20	18x12x7" (46x30x18cm)

Materials

- Non-Woven Green Polypropylene (PP)
- Aluminum Laminate on Polyethylene (PE) Film
- Nylon Hook and Loop

Product Attributes

- High thermal retention
- Non-conductive

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Not Applicable

Product Certifications

- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions

Kit Components


Component	Product Code	Dimensions	Qty
Blanket	68-480151	48x48" (122x122cm)	1
Adult Cap	68-480261	12" (30cm)	1
Leg Warmers	68-480751	16" x 34" (41x86cm)	2



Blanket
68-480151



Adult Cap
68-480261



Leg Warmer
68-480751



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

