

# Temperature Maintenance

## Cozee™ Kit

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
68-480961	Kit, Blanket w/ Adult Cap, 48x84" (122x213cm)	24	18x12x8" (46x30x21cm)

### Materials

- Lamination Fabric of PP with Aluminium Coated on PE Film, Light Green

### Product Attributes


- Three-Layer Warmth – 20% better insulation than cotton blankets
- High thermal retention
- Non-conductive

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.



Product Code	A	B	C	D
68-480961	48" (122cm)	84" (213cm)	20" (51cm)	12" (30cm)

### Product Certifications/Standards

- CE certified as class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

