

Equipment Covers

ECO² Dome Style Band Bags

	CODE	DESCRIPTION (DxWuxWs)	CASE QTY	CASE DIMENSIONS
N	E40-0014	ECO ² Band Bag, Dome Style, Sewn Elastic, 14x9x26" (36x23x66cm)	100	18x12x6" (46x30x14cm)
S	E40-0014-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 14x9x26" (36x23x66cm)	25	13x9x9" (33x23x23cm)
N	E40-0018	ECO ² Band Bag, Dome Style, Sewn Elastic, 18x12x33" (46x30x84cm)	100	18x12x8" (46x30x21cm)
S	E40-0018-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 18x12x33" (46x30x84cm)	25	13x9x13" (32x23x32cm)
N	E40-0020	ECO ² Band Bag, Dome Style, Sewn Elastic, 20x15x37" (51x38x94cm)	100	21x17x7" (54x44x17cm)
S	E40-0020-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 20x15x37" (51x38x94cm)	50	21x17x7" (54x44x17cm)
N	E40-0022	ECO ² Band Bag, Dome Style, Sewn Elastic, 22x15x41" (56x38x104cm)	100	21x17x7" (54x44x17cm)
S	E40-0022-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 22x15x41" (56x38x104cm)	25	13x9x13" (32x23x32cm)
N	E40-0026	ECO ² Band Bag, Dome Style, Sewn Elastic, 26x16x46" (66x41x117cm)	100	15x14x14" (39x36x36cm)
S	E40-0026-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 26x16x46" (66x41x117cm)	25	15x10x14" (37x25x36cm)
N	E40-0027	ECO ² Band Bag, Dome Style, Sewn Elastic, 27x16x46" (69x40x117cm)	100	20x13x14" (52x33x35cm)
S	E40-0027-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 27x16x46" (69x40x117cm)	25	15x10x14" (37x25x36cm)
N	E40-0030	ECO ² Band Bag, Dome Style, Sewn Elastic, 30x16x48" (76x41x122cm)	100	21x13x14" (52x33x35cm)
S	E40-0030-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 30x16x48" (76x41x122cm)	25	15x10x14" (37x25x36cm)
N	E40-0040	ECO ² Band Bag, Dome Style, Sewn Elastic, 40x18x54" (102x46x137cm)	25	16x10x9" (40x25x22cm)
S	E40-0040-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 40x18x54" (102x46x137cm)	25	20x13x16" (51x33x41cm)
N	E40-0050	ECO ² Band Bag, Dome Style, Sewn Elastic, 50x18x54" (127x46x137cm)	25	20x14x9" (50x35x22cm)
S	E40-0050-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 50x18x54" (127x46x137cm)	25	20x13x16" (51x33x41cm)
N	E40-0060	ECO ² Band Bag, Dome Style, Sewn Elastic, 60x16x48" (152x41x122cm)	25	18x11x13" (45x29x33cm)
S	E40-0060-S	ECO ² Band Bag, Dome Style, Sewn Elastic, 60x16x48" (152x41x122cm)	25	18x14x20" (45x35x52cm)

Materials

- ECO² Clear Polyethylene (PE) Film
- Elastic

Product Attributes

- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability
- Dome style shape
- Clear film helps prevent image distortion

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.



Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
 www.alleset.com
 SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.
THINK.MAKE.PROTECT.

GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
 www.alleset.com
 US FDA Establishment Registration #3004911384

EC REP Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
 www.emergobyul.com

