

CORE Room Turnaround Kits®

ECO² Drawcord Bag

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
E50-DCBYL01	ECO ² Drawcord Bag, Yellow, 29x38" (74x97cm)	150	16x12.5x10" (40x32x25cm)
E50-DCBYL02	ECO ² Drawcord Bag, Yellow, 38x38" (97x97cm)	32	TBD

Materials

- Yellow Polyethylene (PE) Film
- Black Polypropylene (PP) Cord

Product Attributes

- Incorporates sustainable materials which enhance product biodegradation
- Drawcord closure with cord lock

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 1 year based upon sterile packaged product aging studies and normal warehouse storage conditions. After 1 year, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage


- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

