

Equipment Covers

Remote Control Cover

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
<input type="radio"/> N 46-RC411	Remote Control Cover, Adhesive Closure, 4x11" (10x28cm)	100	11x9x4" (29x24x11cm)
<input type="radio"/> S 46-RC411-S	Remote Control Cover, Adhesive Closure, 4x11" (10x28cm)	50	11x9x4" (29x24x11cm)
<input type="radio"/> N 46-RC411Z	Remote Control Cover, Zip-tie Closure, 4x11" (10x28cm)	100	11x9x4" (29x24x11cm)
<input type="radio"/> S 46-RC411Z-S	Remote Control Cover, Zip-tie Closure, 4x11" (10x28cm)	50	11x9x4" (29x24x11cm)
<input type="radio"/> N 46-RC512	Remote Control Cover, Adhesive Closure, 5x12" (13x30cm)	100	11x9x4" (29x24x11cm)
<input type="radio"/> S 46-RC512-S	Remote Control Cover, Adhesive Closure, 5x12" (13x30cm)	50	11x9x4" (29x24x11cm)
<input type="radio"/> N 46-RC714	Remote Control Cover, Adhesive Closure, 7x14" (18x36cm)	100	11x9x4" (29x24x11cm)
<input type="radio"/> S 46-RC714-S	Remote Control Cover, Adhesive Closure, 7x14" (18x36cm)	50	11x9x4" (29x24x11cm)
<input type="radio"/> N 46-RC714V	Remote Control Cover, Hook and Loop Closure, 7x14" (18x36cm)	100	11x9x4" (29x24x11cm)
<input type="radio"/> S 46-RC714V-S	Remote Control Cover, Hook and Loop Closure, 7x14" (18x36cm)	50	11x9x4" (29x24x11cm)

Materials


- Clear Polyethylene (PE) Film
- Adhesive Tape, Hook and Loop Material, Zip-tie Material

Product Attributes

- Clear film allows visibility
- Multiple closure configurations

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.



Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleset.com
 SM-TDS-003 2024
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
 US FDA Establishment Registration #3004911384

EC REP
 Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.emergobyul.com

