

# Equipment Covers

## ECO² Remote Control Cover

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
E46-RC411	ECO² Remote Control Cover, Adhesive Closure, 4x11" (10x28cm)	100	11x9x4" (29x24x11cm)
E46-RC411-S	ECO² Remote Control Cover, Adhesive Closure, 4x11" (10x28cm)	50	11x9x4" (29x24x11cm)
E46-RC411Z	ECO² Remote Control Cover, Zip-tie Closure, 4x11" (10x28cm)	100	11x9x4" (29x24x11cm)
E46-RC411Z-S	ECO² Remote Control Cover, Zip-tie Closure, 4x11" (10x28cm)	50	11x9x4" (29x24x11cm)
E46-RC512	ECO² Remote Control Cover, Adhesive Closure, 5x12" (13x30cm)	100	11x9x4" (29x24x11cm)
E46-RC512-S	ECO² Remote Control Cover, Adhesive Closure, 5x12" (13x30cm)	50	11x9x4" (29x24x11cm)
E46-RC714	ECO² Remote Control Cover, Adhesive Closure, 7x14" (18x36cm)	100	11x9x4" (29x24x11cm)
E46-RC714-S	ECO² Remote Control Cover, Adhesive Closure, 7x14" (18x36cm)	50	11x9x4" (29x24x11cm)
E46-RC714V	ECO² Remote Control Cover, Hook and Loop Closure, 7x14" (18x36cm)	100	11x9x4" (29x24x11cm)
E46-RC714V-S	ECO² Remote Control Cover, Hook and Loop Closure, 7x14" (18x36cm)	50	11x9x4" (29x24x11cm)

### Materials

- ECO² Clear Polyethylene (PE) Film
- Adhesive Tape, Hook and Loop Material, Zip-tie Material

### Product Attributes

- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability
- Clear film allows visibility
- Multiple closure configurations

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

### Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
[www.alleset.com](http://www.alleset.com)  
 SM-TDS-003 2023  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
 Westervoortsedijk 60,  
 6827 AT Arnhem,  
 The Netherlands  
[www.emergobyul.com](http://www.emergobyul.com)