


# Equipment Covers

## Dome Style Band Bags



CODE	DESCRIPTION (DxWuxWs)	CASE QTY	CASE DIMENSIONS
<b>N</b> 40-0009	Band Bag, Dome Style, Sewn Elastic, 9x8x14" (23x20x36cm)	100	18x12x6" (46x30x14cm)
<b>S</b> 40-0009-S	Band Bag, Dome Style, Sewn Elastic, 9x8x14" (23x20x36cm)	25	11x9x6" (29x24x16cm)
<b>N</b> 40-0014	Band Bag, Dome Style, Sewn Elastic, 14x9x26" (36x23x66cm)	100	18x12x6" (46x30x14cm)
<b>S</b> 40-0014-S	Band Bag, Dome Style, Sewn Elastic, 14x9x26" (36x23x66cm)	25	13x9x9" (33x23x23cm)
<b>N</b> 40-0018	Band Bag, Dome Style, Sewn Elastic, 18x12x33" (46x30x84cm)	100	18x12x8" (46x30x21cm)
<b>S</b> 40-0018-S	Band Bag, Dome Style, Sewn Elastic, 18x12x33" (46x30x84cm)	25	13x9x13" (32x23x32cm)
<b>N</b> 40-0020	Band Bag, Dome Style, Sewn Elastic, 20x15x37" (51x38x94cm)	100	21x17x7" (54x44x17cm)
<b>S</b> 40-0020-S	Band Bag, Dome Style, Sewn Elastic, 20x15x37" (51x38x94cm)	50	21x17x7" (54x44x17cm)
<b>N</b> 40-0022	Band Bag, Dome Style, Sewn Elastic, 22x15x41" (56x38x104cm)	100	21x17x7" (54x44x17cm)
<b>S</b> 40-0022-S	Band Bag, Dome Style, Sewn Elastic, 22x15x41" (56x38x104cm)	25	13x9x13" (32x23x32cm)
<b>N</b> 40-0026	Band Bag, Dome Style, Sewn Elastic, 26x16x46" (66x41x117cm)	100	15x14x14" (39x36x36cm)
<b>S</b> 40-0026-S	Band Bag, Dome Style, Sewn Elastic, 26x16x46" (66x41x117cm)	25	15x10x14" (37x25x36cm)
<b>N</b> 40-0027	Band Bag, Dome Style, Sewn Elastic, 27x16x46" (69x40x117cm)	100	20x13x14" (52x33x35cm)
<b>S</b> 40-0027-S	Band Bag, Dome Style, Sewn Elastic, 27x16x46" (69x40x117cm)	25	15x10x14" (37x25x36cm)
<b>N</b> 40-0030	Band Bag, Dome Style, Sewn Elastic, 30x16x48" (76x41x122cm)	100	21x13x14" (52x33x35cm)
<b>S</b> 40-0030-S	Band Bag, Dome Style, Sewn Elastic, 30x16x48" (76x41x122cm)	25	15x10x14" (37x25x36cm)
<b>N</b> 40-0040	Band Bag, Dome Style, Sewn Elastic, 40x18x54" (102x46x137cm)	25	16x10x9" (40x25x22cm)
<b>S</b> 40-0040-S	Band Bag, Dome Style, Sewn Elastic, 40x18x54" (102x46x137cm)	25	20x13x16" (51x33x41cm)
<b>N</b> 40-0050	Band Bag, Dome Style, Sewn Elastic, 50x18x54" (127x46x137cm)	25	20x14x9" (50x35x22cm)
<b>S</b> 40-0050-S	Band Bag, Dome Style, Sewn Elastic, 50x18x54" (127x46x137cm)	25	20x13x16" (51x33x41cm)
<b>N</b> 40-0060	Band Bag, Dome Style, Sewn Elastic, 60x16x48" (152x41x122cm)	25	18x11x13" (45x29x33cm)
<b>S</b> 40-0060-S	Band Bag, Dome Style, Sewn Elastic, 60x16x48" (152x41x122cm)	25	18x14x20" (45x35x52cm)

### Materials

- Clear Polyethylene (PE) Film
- Elastic

### Product Attributes

- Dome style shape
- Clear film helps prevent image distortion

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years

**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.  
**THINK.MAKE.PROTECT.**



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

