

Equipment Covers

Rectangular Band Bags

CODE	DESCRIPTION (DxWuxW)	CASE QTY	CASE DIMENSIONS
N 40-2036	Band Bag w/ Tape & Rubber Band, 20x10x36" (51x25x91cm)	100	13x19x13" (33x23x32cm)
S 40-2036-S	Band Bag w/ Tape & Rubber Band, 20x10x36" (51x25x91cm)	25	13x19x13" (33x23x32cm)
N 40-2812	Band Bag w/ Tape & Rubber Band, 28x10x12" (71x25x30cm)	100	18x12x6" (46x30x14cm)
S 40-2812-S	Band Bag w/ Tape & Rubber Band, 28x10x12" (71x25x30cm)	25	13x19x13" (33x23x32cm)
N 40-3020	Band Bag w/ Tape & Rubber Band, 20x10x30" (51x25x76cm)	100	18x12x7" (46x30x18cm)
S 40-3020-S	Band Bag w/ Tape & Rubber Band, 20x10x30" (51x25x76cm)	25	15x10x14" (37x25x36cm)
N 40-3030	Band Bag w/ Tape & Rubber Band, 30x10x30" (76x25x76cm)	100	19x11x13" (47x29x33cm)
S 40-3030-S	Band Bag w/ Tape & Rubber Band, 30x10x30" (76x25x76cm)	25	15x10x14" (37x25x36cm)
N 40-3036	Band Bag w/ Tape & Rubber Band, 36x10x30" (91x25x76cm)	100	19x11x13" (47x29x33cm)
S 40-3036-S	Band Bag w/ Tape & Rubber Band, 36x10x30" (91x25x76cm)	25	15x10x14" (37x25x36cm)
N 40-3620	Band Bag w/ Tape & Rubber Band, 20x10x36" (51x25x91cm)	100	19x11x13" (47x29x33cm)
S 40-3620-S	Band Bag w/ Tape & Rubber Band, 20x10x36" (51x25x91cm)	25	13x9x13" (32x23x32cm)
N 40-3628	Band Bag w/ Tape & Rubber Band, 28x10x36" (71x25x91cm)	100	19x11x13" (47x29x33cm)
S 40-3628-S	Band Bag w/ Tape & Rubber Band, 28x10x36" (71x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3630	Band Bag w/ Tape & Rubber Band, 30x10x36" (76x25x91cm)	100	19x11x13" (47x29x33cm)
S 40-3630-S	Band Bag w/ Tape & Rubber Band, 30x10x36" (76x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3636	Band Bag w/ Tape & Rubber Band, 36x10x36" (91x25x91cm)	100	21x17x7" (54x44x17cm)
S 40-3636-S	Band Bag w/ Tape & Rubber Band, 36x10x36" (91x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3640	Band Bag w/ Tape & Rubber Band, 40x10x36" (102x25x91cm)	100	21x17x7" (54x44x17cm)
S 40-3640-S	Band Bag w/ Tape & Rubber Band, 40x10x36" (102x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3650	Band Bag w/ Tape & Rubber Band, 50x10x36" (127x25x91cm)	100	21x17x8" (54x44x20cm)
S 40-3650-S	Band Bag w/ Tape & Rubber Band, 50x10x36" (127x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3650R	Band Bag w/ Tape & Two Rubber Bands, 50x10x36" (127x25x91cm)	100	21x17x8" (54x44x20cm)
S 40-3650R-S	Band Bag w/ Tape & Two Rubber Bands, 50x10x36" (127x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3654	Band Bag w/ Tape & Rubber Band, 54x10x36" (137x25x91cm)	100	21x17x8" (54x44x20cm)
S 40-3654-S	Band Bag w/ Tape & Rubber Band, 54x10x36" (137x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-3654R	Band Bag w/ Tape & Two Rubber Bands, 54x10x36" (137x25x91cm)	100	21x17x8" (54x44x20cm)
S 40-3654R-S	Band Bag w/ Tape & Two Rubber Bands, 54x10x36" (137x25x91cm)	25	15x10x14" (37x25x36cm)
N 40-4020	Band Bag w/ Tape, Sewn Elastic, 20x14x40" (51x36x102cm)	100	21x17x7" (54x44x17cm)
S 40-4020-S	Band Bag w/ Tape, Sewn Elastic, 20x14x40" (51x36x102cm)	25	15x10x14" (37x25x36cm)
N 40-4030	Band Bag w/ Tape, Sewn Elastic, 30x14x40" (76x36x102cm)	100	21x17x7" (54x44x17cm)
S 40-4030-S	Band Bag w/ Tape, Sewn Elastic, 30x14x40" (76x36x102cm)	25	15x10x14" (37x25x36cm)
N 40-4040	Band Bag w/ Tape, Sewn Elastic, 40x14x40" (102x36x102cm)	100	21x17x8" (54x44x20cm)
S 40-4040-S	Band Bag w/ Tape, Sewn Elastic, 40x14x40" (102x36x102cm)	25	15x10x14" (37x25x36cm)
N 40-4050	Band Bag w/ Tape, Sewn Elastic, 30x14x40" (76x36x102cm)	100	21x17x8" (54x44x20cm)
N 40-4436	Band Bag w/ Tape, Sewn Elastic, 36x15x44" (91x38 x112cm)	100	21x17x6" (54x44x14cm)
S 40-4436-S	Band Bag w/ Tape, Sewn Elastic, 36x15x44" (91x38 x112cm)	25	15x10x14" (37x25x36cm)
N 40-5030	Band Bag w/ Tape, Sewn Elastic, 30x17x50" (76x43x127cm)	100	21x17x8" (54x44x20cm)
S 40-5030-S	Band Bag w/ Tape, Sewn Elastic, 30x17x50" (76x43x127cm)	25	15x10x14" (37x25x36cm)
N 40-5040	Band Bag w/ Tape, Sewn Elastic, 40x17x50" (102x43x127cm)	100	21x17x8" (54x44x20cm)
S 40-5040-S	Band Bag w/ Tape, Sewn Elastic, 40x17x50" (102x43x127cm)	25	15x10x14" (37x25x36cm)
S 40-3028-S	Band Bag w/ Tape & Rubber Band, 28x10x30" (71x25x76cm)	25	15x10x14" (37x25x36cm)
N 40-3040	Band Bag w/ Tape & Rubber Band, 40x10x30" (102x25x76cm)	100	19x11x13" (47x29x33cm)
N 40-5050	Band Bag w/ Tape & Rubber Band, 50x10x50" (127x25x127cm)	100	21x17x8" (54x44x20cm)
S 40-5050-S	Band Bag w/ Tape & Rubber Band, 50x10x50" (127x25x127cm)	25	15x10x14" (37x25x36cm)
N 40-2020	Band Bag w/ Tape & Rubber Band, 20x10x20" (51x25x51cm)	100	11x9x10" (28x22x25cm)
S 40-2020-S	Band Bag w/ Tape & Rubber Band, 20x10x20" (51x25x51cm)	25	21x17x7" (54x44x17cm)
N 40-2028	Band Bag w/ Tape & Rubber Band, 20x10x28" (51x25x71cm)	100	11x9x10" (28x22x25cm)
S 40-2028-S	Band Bag w/ Tape & Rubber Band, 20x10x28" (51x25x71cm)	25	21x17x7" (54x44x17cm)

Equipment Covers

Rectangular Band Bags

CODE	DESCRIPTION (DxWuxW)	CASE QTY	CASE DIMENSIONS
N 40-2820	Band Bag w/ Tape & Rubber Band, 28x10x20" (51x25x51cm)	100	21x17x7" (54x44x17cm)
S 40-2820-S	Band Bag w/ Tape & Rubber Band, 28x10x20" (51x25x51cm)	25	13x9x13" (33x23x32cm)
N 40-2836	Band Bag w/ Tape & Rubber Band, 28x10x36" (71x25x91cm)	100	21x17x7" (54x44x17cm)
S 40-2836-S	Band Bag w/ Tape & Rubber Band, 28x10x36" (71x25x91cm)	25	13x9x13" (33x23x32cm)
N 40-3028	Band Bag w/ Tape & Rubber Band, 30x10x28" (76x25x71cm)	100	21x17x7" (54x44x17cm)
S 40-3040-S	Band Bag w/ Tape & Rubber Band, 30x10x40" (76x25x102cm)	25	13x9x13" (33x23x32cm)
N 40-3612	Band Bag w/ Tape & Rubber Band, 36x10x12" (91x25x30cm)	100	21x17x7" (54x44x17cm)
S 40-3612-S	Band Bag w/ Tape & Rubber Band, 36x10x12" (91x25x30cm)	25	13x9x13" (33x23x32cm)
N 40-1212	Band Bag w/ Tape & Rubber Band, 12x10x12" (30x25x30cm)	100	19x11x13" (47x29x33cm)

Materials

- Clear Polyethylene (PE) Film
- Elastic & Rubber Band

Product Attributes

- Rectangular shape
- Clear film helps prevent image distortion
- Elastic closure
- Rubber band securement (select covers)

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleset.com
 SM-TDS-003 2024
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.emergobyul.com