

Surgical Drapes

Mayo Stand Cover

CODE	DESCRIPTION (WxL)	CASE QTY	CASE DIMENSIONS
80-13105	Mayo Stand Cover, 23x54" (58x137cm)	180	22x17x9" (55x44x23cm)
80-13105-S	Mayo Stand Cover, 23x54" (58x137cm)	60	19x13x12" (49x33x30cm)
80-13106	Mayo Stand Cover, 30x57" (76x145cm)	140	18x12x19" (46x30x49cm)
80-13106-S	Mayo Stand Cover, 30x57" (76x145cm)	60	21x13x15" (53x33x37cm)

Materials

- Multilayered SMS Material with Antistatic Properties
- Blue Polyethylene (PE) Film

Product Attributes

- SMS reinforced fabric

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleaset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleaset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

