

Equipment Covers

ECO² X-Ray Cassette

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
E46-CAS2024L	ECO ² X-Ray Cassette Cover, Medium, Cuffed Release Liner, 20x24" (51x61cm)	100	20x11x19" (52x29x47cm)
E46-CAS2024L-S	ECO ² X-Ray Cassette Cover, Medium, Cuffed Release Liner, 20x24" (51x61cm)	50	20x11x19" (52x29x47cm)

Materials

- ECO² Clear Polyethylene (PE) Film
- ECO² Blue Polyethylene (PE) Film

Product Attributes

- Clear Polyethylene film allows for visibility
- Cuffed opening
- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

