

Equipment Covers

Foot Switch Covers

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
 46-FS1424	Foot Switch Cover, 14x24" (35.5x61cm)	100	18x12x7" (46x30x18cm)
 46-FS1430	Foot Switch Cover, 14x30" (35.5x61cm)	100	18x12x7" (46x30x18cm)
 46-FS1618	Foot Switch Cover, 16x18" (41x46cm)	100	14x11x8" (36x29x21cm)
 46-FS2020	Foot Switch Cover, 20x20" (51x51cm)	100	18x12x8" (46x30x21cm)
 46-FS2420	Foot Switch Cover, 24x20" (61x51cm)	100	18x12x8" (46x30x21cm)
 46-FS3030	Foot Switch Cover, 30x30" (76x76cm)	100	18x16x10" (45x41x26cm)

Materials

- Clear Polyethylene (PE) Film
- Polypropylene (PP) Filament String

Product Attributes

- Drawstring closure
- Clear film allows for visibility

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified


Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

 GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

 Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com