

Heel Incision

gentleheel®

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
MPGH10X50	gentleheel® Micro-Premie, Yellow, 1.5x1.3" (3.9x3.2cm)	50EA/Pack, 10Pack/CS	11x11x17" (28x28x43cm)
NGH10X50	gentleheel® Newborn, Green, 1.5x1.3" (3.9x3.2cm)	50EA/Pack, 10Pack/CS	11x11x17" (28x28x43cm)
PGH10X50	gentleheel® Premie, Lavender, 1.5x1.3" (3.9x3.2cm)	50EA/Pack, 10Pack/CS	11x11x17" (28x28x43cm)
TGH10X50	gentleheel® Toddler, Blue, 1.5x1.3" (3.9x3.2cm)	50EA/Pack, 10Pack/CS	11x11x17" (28x28x43cm)

Materials

- Housing/Button/Locker: Acrylonitrile Butadiene Styrene (ABS)
- Blade holder/Spring holder: Polyetherimide
- Blade: Stainless Steel
- Spring: Nickel plating spring steels

Product Attributes

- Color coded for easy identification
- Designed to be easy to use
- Incision point clearly marked

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Shelf Life

Sterile: 5 years

Product Certifications/Standards


- CE certified as Class IIa, according to MDR 2017/745 Annex VIII, Rule 6

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage


- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.alleset.com

