

# Gowns

## Non-Reinforced

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
<b>N</b> V90-2010	Gown, Non-Reinforced	M	48	18x15x12" (46x38x29.5cm)
<b>S</b> V90-2010-S2	Gown, Non-Reinforced, 2 Hand Towels	M	28	16x11x18" (42x28x45cm)
<b>N</b> V90-2020	Gown, Non-Reinforced	L	48	19x16x12.5" (50x42x32cm)
<b>S</b> V90-2020-S2	Gown, Non-Reinforced, 2 Hand Towels	L	26	16x11x18" (42x28x45cm)
<b>N</b> V90-2021	Gown, Non-Reinforced	L, X-Long	48	19x16x12.5" (50x42x32cm)
<b>S</b> V90-2021-S2	Gown, Non-Reinforced, 2 Hand Towels	L, X-Long	26	16x11x18" (42x28x45cm)
<b>N</b> V90-2030	Gown, Non-Reinforced	XL	48	19x16x12.5" (50x42x32cm)
<b>S</b> V90-2030-S2	Gown, Non-Reinforced, 2 Hand Towels	XL	24	16x11x18" (42x28x45cm)
<b>N</b> V90-2031	Gown, Non-Reinforced	XL, X-Long	48	19x16x12.5" (50x42x32cm)
<b>S</b> V90-2031-S2	Gown, Non-Reinforced, 2 Hand Towels	XL, X-Long	24	16x11x18" (42x28x45cm)
<b>N</b> V90-2040	Gown, Non-Reinforced	2XL	44	19x16x12.5" (50x42x32cm)
<b>S</b> V90-2040-S2	Gown, Non-Reinforced, 2 Hand Towels	2XL	24	16x11x18" (42x28x45cm)

### Materials

- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency

### Product Attributes

- Sonically bonded seams throughout gown
- Adjustable hook and loop fenestration
- Knit cuffs
- Set-in sleeves

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- ANSI/AAMI PB70:2012 Level 2 Liquid Barrier Protection
- EN 13795:2011 Standard Performance
- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 4
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Code	A: Back Length	B: Sleeve Length
V90-2010 (S2/NS)	41" (104cm)	23" (58cm)
V90-2020 (S2/NS)	45" (115cm)	24" (62cm)
V90-2021(S2/NS)	54" (138cm)	24" (62cm)
V90-2030 (S2/NS)	49" (124cm)	25" (64cm)
V90-2031(S2/NS)	58" (147cm)	25" (64cm)
V90-2040 (S2/NS)	51" (130cm)	30" (75cm)
Tolerance	±1" (2.54cm)	±0.5" (1.3cm)



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
[www.alleset.com](http://www.alleset.com)  
 SM-TDS-002 2022  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
 US FDA Establishment Registration #3004911384

**EC REP**  
 Emergo Europe B.V.  
 Prinsessegracht 20  
 2514 AP The Hague  
 The Netherlands  
[www.emergobyul.com](http://www.emergobyul.com)