

# Patient Positioning

## Head Cradles

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
50-150101	Head Cradle, Adult, 9x8x4.25" (23x20x10.8cm)	36	24x17x25" (61x43x64cm)
50-150102	Head Cradle with Trach Slit, Adult, 9x8x4.25" (23x20x10.8cm)	36	24x17x25" (61x43x64cm)
50-150104	Head Cradle, Stacked with Trach Slit, Adult, 8x8.5x4" (20x22x10cm)	36	24x17x25" (61x43x64cm)
50-150105	Head Rest, Donut, 9x2" (23x5cm)	36	20.5x18x19" (52x46x49cm)
50-150106	Head Rest, Donut Multiple, 9x2" (23x5cm)	36	20.5x18x19" (52x46x49cm)
50-150109	Head Cradle, Pediatric, 6.25x7x3.5" (16x18x8cm)	36	20.5x18x19" (52x46x49cm)

### Materials

- Polyurethane Foam

### Product Attributes

- High density polyurethane foam

### Sterilization

Non-Sterile: This product is intended to be used non-sterile. If the repackager intends to sterilize this product, a validation must be conducted to demonstrate their sterilization process will meet all the ISO 11135 and EN 10993 requirements to demonstrate sterility and EO/ECH residual compliance. It must be done according to their validated sterilization process.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Alleaset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
 www.alleaset.com  
 SM-TDS-002 2022  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

Global Resources International Tennessee Foam Inc,  
 910 Callahan Dr #108,  
 Knoxville, Tennessee,  
 37912 UNITED STATES  
 www.alleaset.com  
 US FDA Establishment Registration #3017728970

Alleaset Limited B.V.  
 Kingsfordweg 151,  
 1043 GR Amsterdam,  
 The Netherlands  
 www.alleaset.com