

Surgical Accessories

Light Handle Covers

| CODE | DESCRIPTION | CASE QTY | CASE DIMENSIONS |
|------------|-----------------------------|----------|------------------------|
| 25-3011 | Flexible Light Handle Cover | 2400 | 31x24x20" (79x61x51cm) |
| 25-3011-S1 | Flexible Light Handle Cover | 128 | 15x12x8" (38x30x20cm) |
| 25-3011-S2 | Flexible Light Handle Cover | 160 | 15x12x8" (38x30x20cm) |
| 25-3011-S3 | Flexible Light Handle Cover | 240 | 15x12x8" (38x30x20cm) |

Materials

- Low Density Polyethylene (LDPE)


Product Attributes

- Flexible, fits most surgical light handle adapters
- Reduces need to reprocess and sterilize

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.



Shelf Life

Sterile: 5 Years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, product should be performance tested to ensure compliance to specifications

Product Certifications/Standards

- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands.
www.emergobyul.com

