

CORE Room Turnaround Kits®

ECO² Liner/Miscellaneous Bag

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
E50-CTB01	ECO ² Bag, Clear, 14x18" (36x46cm)	TBD	TBD
E50-CTB02	ECO ² Bag, Clear, 14x24" (36x61cm)	TBD	TBD
E50-CTB03	ECO ² Bag, Clear, 16x9x30" (41x23x76cm)	250	12x9x10" (31x23x25cm)
E50-CTB05	ECO ² Bag, Clear, 38x60" (97x152cm)	150	16x11x11" (40x28x29cm)
E50-CTB06	ECO ² Trash Liner Bag, Clear, 33x38 " (84x97cm)	250	16x12.5x10" (40x32x25cm)
E50-CTB07	ECO ² Bag, Red, 21x19x46" (53x48x117cm)	TBD	TBD
E50-CTB08	ECO ² Bag, Clear, 40x46" (102x117cm)	TBD	TBD
E50-CTB09	ECO ² Bag, Brown, 38x38" (97x97cm)	100	16x12.5x10" (40x32x25cm)
E50-CTB11	ECO ² Bag, Brown 38x58" (97x147cm)	TBD	TBD

Materials

Polyethylene (PE) Film

Product Attributes

- Available in various film colors and sizes
- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 1 year based upon sterile packaged product aging studies and normal warehouse storage conditions. After 1 year, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards



- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2024
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

GRI-Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.emergobyul.com

