

Patient Positioning

Trendelenburg Position System

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
50-580	Trendelenburg Positioning System, 29x20" (74x51cm)	6 set/case	31x22x9" (80x55x24cm)
50-582	Trendelenburg Positioning System, 35x20" (89x51cm)	6 set/case	37x22x9" (95x55x24cm)
50-583	Trendelenburg Positioning System, 40x20" (102x51cm)	6 set/case	43x22x10" (108x55x26cm)

Materials

- Table Pad: Polyurethane (PU) Foam, Green
- Straps: Composite Sponge and Velcro Hook, White
- Lift Sheet: Polyethylene Terephthalate (PET) Spunlace

Product Attributes

- Viscoelastic "memory" foam pad
- Integrated pad securement straps
- Foam patient positioning straps
- Spunlace lift sheet for patient transfer

Sterilization

Non-Sterile: This product is intended to be used non-sterile. If the repackager intends to sterilize this product, a validation must be conducted to demonstrate their sterilization process will meet all the ISO 11135 and EN 10993 requirements to demonstrate sterility and EO/ECH residual compliance. It must be done according to their validated sterilization process.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-0022022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

