


Equipment Covers

ECO² Rectangular Band Bags

CODE	DESCRIPTION (DxW)	CASE QTY	CASE DIMENSIONS
N E40-1212	ECO ² Band Bag w/ Tape & Rubber Band, 12x12" (30x30cm)	100	18x12x5" (46x30x14cm)
N E40-2020	ECO ² Band Bag w/ Tape & Rubber Band, 20x20" (51x51cm)	100	21x17x7" (54x44x17cm)
S E40-2020-S	ECO ² Band Bag w/ Tape & Rubber Band, 20x20" (51x51cm)	25	11x9x10" (28x22x25cm)
N E40-2028	ECO ² Band Bag w/ Tape & Rubber Band, 20x28" (51x71cm)	100	21x17x7" (54x44x17cm)
S E40-2028-S	ECO ² Band Bag w/ Tape & Rubber Band, 20x28" (51x71cm)	25	11x9x10" (28x22x25cm)
N E40-2036	ECO ² Band Bag w/ Tape & Rubber Band, 20x36" (51x91cm)	100	13x19x13" (33x23x32cm)
S E40-2036-S	ECO ² Band Bag w/ Tape & Rubber Band, 20x36" (51x91cm)	25	13x19x13" (33x23x32cm)
N E40-2812	ECO ² Band Bag w/ Tape & Rubber Band, 28x12" (71x30cm)	100	18x12x6" (46x30x14cm)
S E40-2812-S	ECO ² Band Bag w/ Tape & Rubber Band, 28x12" (71x30cm)	25	13x19x13" (33x23x32cm)
N E40-2820	ECO ² Band Bag w/ Tape & Rubber Band, 28x20" (71x51cm)	100	21x17x7" (54x44x17cm)
S E40-2820-S	ECO ² Band Bag w/ Tape & Rubber Band, 28x20" (71x51cm)	25	13x9x13" (33x23x32cm)
N E40-2836	ECO ² Band Bag w/ Tape & Rubber Band, 28x36" (71x91cm)	100	19x11x13" (47x29x33cm)
S E40-2836-S	ECO ² Band Bag w/ Tape & Rubber Band, 28x36" (71x91cm)	25	15x10x14" (37x25x36cm)
N E40-3020	ECO ² Band Bag w/ Tape & Rubber Band, 30x20" (76x51cm)	100	18x12x7" (46x30x18cm)
S E40-3020-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x20" (76x51cm)	25	15x10x14" (37x25x36cm)
N E40-3028	ECO ² Band Bag w/ Tape & Rubber Band, 30x28" (76x71cm)	100	19x11x13" (47x29x33cm)
S E40-3028-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x28" (76x71cm)	25	15x10x14" (37x25x36cm)
N E40-3030	ECO ² Band Bag w/ Tape & Rubber Band, 30x30" (76x76cm)	100	19x11x13" (47x29x33cm)
S E40-3030-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x30" (76x76cm)	25	15x10x14" (37x25x36cm)
N E40-3036	ECO ² Band Bag w/ Tape & Rubber Band, 30x36" (76x91cm)	100	19x11x13" (47x29x33cm)
S E40-3036-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x36" (76x91cm)	25	15x10x14" (37x25x36cm)
N E40-30E40	ECO ² Band Bag w/ Tape & Rubber Band, 30x28" (76x71cm)	100	19x11x13" (47x29x33cm)
S E40-30E40-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x40" (76x102cm)	25	15x10x14" (37x25x36cm)
N E40-3612	ECO ² Band Bag w/ Tape & Rubber Band, 36x12" (91x30cm)	100	19x11x13" (47x29x33cm)
S E40-3612-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x12" (91x30cm)	25	13x9x13" (33x23x32cm)
N E40-3620	ECO ² Band Bag w/ Tape & Rubber Band, 36x20" (91x51cm)	100	19x11x13" (47x29x33cm)
S E40-3620-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x20" (91x51cm)	25	13x9x13" (32x23x32cm)
N E40-3628	ECO ² Band Bag w/ Tape & Rubber Band, 36x28" (91x71cm)	100	19x11x13" (47x29x33cm)
S E40-3628-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x28" (91x71cm)	25	15x10x14" (37x25x36cm)
N E40-3630	ECO ² Band Bag w/ Tape & Rubber Band, 36x30" (91x76cm)	100	19x11x13" (47x29x33cm)
S E40-3630-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x30" (91x76cm)	25	15x10x14" (37x25x36cm)
N E40-3636	ECO ² Band Bag w/ Tape & Rubber Band, 36x36" (91x91cm)	100	21x17x7" (54x44x17cm)
S E40-3636-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x36" (91x91cm)	25	15x10x14" (37x25x36cm)
N E40-36E40	ECO ² Band Bag w/ Tape & Rubber Band, 36x40" (91x102cm)	100	21x17x7" (54x44x17cm)
S E40-36E40-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x40" (91x102cm)	25	15x10x14" (37x25x36cm)
N E40-3650	ECO ² Band Bag w/ Tape & Rubber Band, 36x50" (91x127cm)	100	21x17x8" (54x44x20cm)
S E40-3650-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x50" (91x127cm)	25	15x10x14" (37x25x36cm)
N E40-3650R	ECO ² Band Bag w/ Tape & Rubber Band, 36x50" (91x127cm)	100	21x17x8" (54x44x20cm)
S E40-3650R-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x50" (91x127cm)	25	15x10x14" (37x25x36cm)
N E40-3654	ECO ² Band Bag w/ Tape & Rubber Band, 36x54" (91x137cm)	100	21x17x8" (54x44x20cm)
S E40-3654-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x54" (91x137cm)	25	15x10x14" (37x25x36cm)
N E40-3654R	ECO ² Band Bag w/ Tape & Rubber Band, 36x54" (91x137cm)	100	21x17x8" (54x44x20cm)
S E40-3654R-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x54" (91x137cm)	25	15x10x14" (37x25x36cm)
N E40-E4020	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x20" (102x51cm)	100	21x17x7" (54x44x17cm)
S E40-E4020-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x20" (102x51cm)	25	15x10x14" (37x25x36cm)
N E40-E4030	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x30" (102x76cm)	100	21x17x7" (54x44x17cm)
S E40-E4030-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x30" (102x76cm)	25	15x10x14" (37x25x36cm)





Equipment Covers

ECO² Rectangular Band Bags

CODE	DESCRIPTION (DxW)	CASE QTY	CASE DIMENSIONS
N E40-4040	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x40" (102x102cm)	100	19x11x13" (47x29x33cm)
S E40-4040-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x40" (102x102cm)	25	15x10x14" (37x25x36cm)
N E40-4050	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x50" (102x127cm)	100	19x11x13" (47x29x33cm)
S E40-4050-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x50" (102x127cm)	25	13x19x13" (33x23x32cm)
N E40-4436	ECO ² Band Bag w/ Tape, Sewn Elastic, 44x36" (112x91cm)	100	21x17x6" (54x44x14cm)
S E40-4436-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 44x36" (112x91cm)	25	15x10x14" (37x25x36cm)
N E40-5030	ECO ² Band Bag w/ Tape, Sewn Elastic, 50x30" (127x76cm)	100	21x17x8" (54x44x20cm)
S E40-5030-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 50x30" (127x76cm)	25	15x10x14" (37x25x36cm)
N E40-5040	ECO ² Band Bag w/ Tape, Sewn Elastic, 50x40" (127x102cm)	100	19x11x13" (47x29x33cm)
S E40-50E40-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 50x40" (127x102cm)	25	15x10x14" (37x25x36cm)
N E40-5050	ECO ² Band Bag w/ Tape, Sewn Elastic, 50x50" (127x127cm)	100	19x11x13" (47x29x33cm)
S E40-5050-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 50x50" (127x127cm)	25	15x10x14" (37x25x36cm)

Materials

- ECO² Clear Polyethylene (PE) Film
- Elastic & Rubber Band

Product Attributes

- Rectangular shape
- Clear film helps prevent image distortion
- Elastic Closure
- Rubber band securement (select covers)
- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability.

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleteset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleteset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleteset.com
 SM-TDS-003 2024
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaying, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleteset.com
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.emergobyul.com

