

Equipment Covers

ECO2 Rectangular Band Bags

	CODE	DESCRIPTION (DxWuxW)	CASE QTY	CASE DIMENSIONS
<input checked="" type="radio"/>	E40-2812	ECO ² Band Bag w/ Tape & Rubber Band	TBD	TBD
<input checked="" type="radio"/>	E40-2812-S	ECO ² Band Bag w/ Tape & Rubber Band	TBD	TBD
<input checked="" type="radio"/>	E40-3020	ECO ² Band Bag w/ Tape & Rubber Band, 20x10x30" (51x25x76cm)	100	18x12x7" (46x30x18cm)
<input checked="" type="radio"/>	E40-3020-S	ECO ² Band Bag w/ Tape & Rubber Band, 20x10x30" (51x25x76cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3030	ECO ² Band Bag w/ Tape & Rubber Band, 30x10x30" (76x25x76cm)	100	19x11x13" (47x29x33cm)
<input checked="" type="radio"/>	E40-3030-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x10x30" (76x25x76cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3036	ECO ² Band Bag w/ Tape & Rubber Band, 36x10x30" (91x25x76cm)	100	19x11x13" (47x29x33cm)
<input checked="" type="radio"/>	E40-3036-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x10x30" (91x25x76cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3620	ECO ² Band Bag w/ Tape & Rubber Band, 20x10x36" (51x25x91cm)	100	19x11x13" (47x29x33cm)
<input checked="" type="radio"/>	E40-3620-S	ECO ² Band Bag w/ Tape & Rubber Band, 20x10x36" (51x25x91cm)	25	13x9x13" (32x23x32cm)
<input checked="" type="radio"/>	E40-3628	ECO ² Band Bag w/ Tape & Rubber Band, 28x10x36" (71x25x91cm)	100	19x11x13" (47x29x33cm)
<input checked="" type="radio"/>	E40-3628-S	ECO ² Band Bag w/ Tape & Rubber Band, 28x10x36" (71x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3630	ECO ² Band Bag w/ Tape & Rubber Band, 30x10x36" (76x25x91cm)	100	19x11x13" (47x29x33cm)
<input checked="" type="radio"/>	E40-3630-S	ECO ² Band Bag w/ Tape & Rubber Band, 30x10x36" (76x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3636	ECO ² Band Bag w/ Tape & Rubber Band, 36x10x36" (91x25x91cm)	100	21x17x7" (54x44x17cm)
<input checked="" type="radio"/>	E40-3636-S	ECO ² Band Bag w/ Tape & Rubber Band, 36x10x36" (91x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3640	ECO ² Band Bag w/ Tape & Rubber Band, 40x10x36" (102x25x91cm)	100	21x17x7" (54x44x17cm)
<input checked="" type="radio"/>	E40-3640-S	ECO ² Band Bag w/ Tape & Rubber Band, 40x10x36" (102x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3650	ECO ² Band Bag w/ Tape & Rubber Band, 50x10x36" (127x25x91cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-3650-S	ECO ² Band Bag w/ Tape & Rubber Band, 50x10x36" (127x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3650R	ECO ² Band Bag w/ Tape & Two Rubber Bands, 50x10x36" (127x25x91cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-3650R-S	ECO ² Band Bag w/ Tape & Two Rubber Bands, 50x10x36" (127x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3654	ECO ² Band Bag w/ Tape & Rubber Band, 54x10x36" (137x25x91cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-3654-S	ECO ² Band Bag w/ Tape & Rubber Band, 54x10x36" (137x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-3654R	ECO ² Band Bag w/ Tape & Two Rubber Bands, 54x10x36" (137x25x91cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-3654R-S	ECO ² Band Bag w/ Tape & Two Rubber Bands, 54x10x36" (137x25x91cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-4020	ECO ² Band Bag w/ Tape, 20x14x37" (51x36x94cm)	100	21x17x7" (54x44x17cm)
<input checked="" type="radio"/>	E40-4020-S	ECO ² Band Bag w/ Tape, 20x14x37" (51x36x94cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-4030	ECO ² Band Bag w/ Tape, Sewn Elastic, 30x14x37" (76x36x94cm)	100	21x17x7" (54x44x17cm)
<input checked="" type="radio"/>	E40-4030-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 30x14x37" (76x36x94cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-4040	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x14x37" (102x36x94cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-4040-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x14x37" (102x36x94cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-4436	ECO ² Band Bag w/ Tape, Sewn Elastic, 36x15x40" (91x38 x102cm)	100	21x17x6" (54x44x14cm)
<input checked="" type="radio"/>	E40-4436-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 36x15x40" (91x38 x102cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-5030	ECO ² Band Bag w/ Tape, Sewn Elastic, 30x17x46" (76x43x117cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-5030-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 30x17x46" (76x43x117cm)	25	15x10x14" (37x25x36cm)
<input checked="" type="radio"/>	E40-5040	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x17x46" (102x43x117cm)	100	21x17x8" (54x44x20cm)
<input checked="" type="radio"/>	E40-5040-S	ECO ² Band Bag w/ Tape, Sewn Elastic, 40x17x46" (102x43x117cm)	25	15x10x14" (37x25x36cm)

Materials

- ECO² Clear Polyethylene (PE) Film
- Elastic & Rubber Band

Product Attributes

- Rectangular shape
- Clear film helps prevent image distortion
- Elastic Closure
- Rubber band securement (select covers)
- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability.



Equipment Covers

ECO2 Rectangular Band Bags

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleset.com
 SM-TDS-003 2023
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.alleset.com

