

RoyalGuard™

Breathable High Performance Surgical Gowns

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
N 91-6110	Breathable High Performance Surgical Gown	S/M	40	20x17x13" (50x42x32cm)
S 91-6110-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	S/M	30	18x12x22" (46x30x55cm)
S 91-6110-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	S/M	30	18x12x22" (46x30x55cm)
N 91-6120	Breathable High Performance Surgical Gown	L	40	20x17x13" (50x42x32cm)
S 91-6120-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	L	28	18x12x22" (46x30x55cm)
S 91-6120-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	L	28	18x12x22" (46x30x55cm)
N 91-6120-HB	Breathable High Performance Surgical Gown, Handi-Bin	L	405	47x31x31" (119x79x78cm)
N 91-6120-W1-HB	Breathable High Performance Surgical Gown, 1 Hand Towel, Handi-Bin	L	378	47x31x31" (119x79x78cm)
N 91-6122	Breathable High Performance Surgical Gown	L, Long	32	20x17x13" (50x42x32cm)
S 91-6122-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	L, Long	26	18x12x22" (46x30x55cm)
S 91-6122-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	L, Long	26	18x12x22" (46x30x55cm)
N 91-6124	Breathable High Performance Surgical Gown, A-Line	L, Long	32	20x17x13" (50x42x32cm)
S 91-6124-S1	Breathable High Performance Surgical Gown, A-Line, 1 Hand Towel	L, Long	24	18x12x22" (46x30x55cm)
S 91-6124-S2	Breathable High Performance Surgical Gown, A-Line, 2 Hand Towels	L, Long	24	18x12x22" (46x30x55cm)
N 91-6130	Breathable High Performance Surgical Gown	XL	40	20x17x13" (50x42x32cm)
S 91-6130-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	XL	26	18x12x22" (46x30x55cm)
S 91-6130-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	XL	26	18x12x22" (46x30x55cm)
N 91-6130-HB	Breathable High Performance Surgical Gown, Handi-Bin	XL	378	47x31x31" (119x79x78cm)
N 91-6130-W1-HB	Breathable High Performance Surgical Gown, 1 Hand Towel, Handi-Bin	XL	351	47x31x31" (119x79x78cm)
N 91-6132	Breathable High Performance Surgical Gown	XL, Long	32	20x17x13" (50x42x32cm)
S 91-6132-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	XL, Long	24	18x12x22" (46x30x55cm)
S 91-6132-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	XL, Long	24	18x12x22" (46x30x55cm)
N 91-6134	Breathable High Performance Surgical Gown, A-Line	XL, Long	32	20x17x13" (50x42x32cm)
S 91-6134-S1	Breathable High Performance Surgical Gown, A-Line, 1 Hand Towel	XL, Long	22	18x12x22" (46x30x55cm)
S 91-6134-S2	Breathable High Performance Surgical Gown, A-Line, 2 Hand Towels	XL, Long	22	18x12x22" (46x30x55cm)
N 91-6140	Breathable High Performance Surgical Gown	2XL	36	20x17x13" (50x42x32cm)
S 91-6140-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	2XL	24	18x12x22" (46x30x55cm)
S 91-6140-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	2XL	24	18x12x22" (46x30x55cm)
N 91-6140-HB	Breathable High Performance Surgical Gown, Handi-Bin	2XL	351	47x31x31" (119x79x78cm)
N 91-6142	Breathable High Performance Surgical Gown	2XL, Long	32	20x17x13" (50x42x32cm)
S 91-6142-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	2XL, Long	22	18x12x22" (46x30x55cm)
N 91-6150	Breathable High Performance Surgical Gown	3XL	32	20x17x13" (50x42x32cm)
S 91-6150-S1	Breathable High Performance Surgical Gown, 1 Hand Towel	3XL	20	18x12x22" (46x30x55cm)
S 91-6150-S2	Breathable High Performance Surgical Gown, 2 Hand Towels	3XL	20	18x12x22" (46x30x55cm)

RoyalGuard™

Breathable High Performance Surgical Gowns

Materials

- Tri-Layer SFS High Performance Breathable Viral Barrier Material
- Tri-Layer SMS Material (ComfortBack™ Panels)

Product Attributes

- SureSleeve™ Sleeve In Sleeve design for maximum protection
- Color coded neckline (red) identifies protection level
- Designed with sonic seaming to eliminate holes caused by stitching
- Protective breathable material
- Adjustable hook and loop neck closure

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards



- EN13795-1:2019 High Performance
- ANSI/AAMI PB70:2012 Level 4 Liquid Barrier Protection
- Class 1 Flammability Standards per 16 CFR 1610
- I1 Protection against laser ignition per ISO 11810
- According to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



A-line

Code	A: Total Length	B: Sleeve Length
91-6110 (S1/S2/NS)	47" (120cm)	21" (54cm)
91-6120 (S1/S2/NS)	49.7" (126cm)	21.5" (54.6cm)
91-6122 (S1/S2/NS)	56" (142cm)	21.5" (54.6cm)
91-6124 (S1/S2/NS)	56.7" (144cm)	22" (56cm)
91-6130 (S1/S2/NS)	52" (132cm)	24.6" (62.5cm)
91-6132 (S1/S2/NS)	58" (148cm)	24.6" (62.5cm)
91-6134 (S1/S2/NS)	60.6" (154cm)	24" (62cm)
91-6140(S1/S2/NS)	52.7" (134cm)	25" (64cm)
91-6142 (S1/NS)	59" (150cm)	25" (64cm)
91-6150 (S1/S2/NS)	56.7" (144cm)	25.6 (65cm)
Tolerance	±0.8" (2cm)	±0.4" (1cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

